


01-8090  
01-8095



**SWL: 200 kg/440 lbs**

## Product Ordering Information

---

Order Ref	Product	Size	
01-8090	FootStool	32 x 38 x 10 cm /12,5 x 15 x 4 inch	
01-8095	FootStool	32 x 38 x 5 cm /12,5 x 15 x 2 inch	



[info@amilakesouthern.co.uk](mailto:info@amilakesouthern.co.uk)

FootStool is used to raise the position of the knees and thereby enable safer, easier bed-to-wheelchair transfer in a seated position. FootStool also provides good support under the feet during visits to the toilet, or when repositioning farther back in a chair. FootStool is made of a material with closed cells, which makes it easy to keep clean.

SystemRoMedic is a series of different models of transfer assistive devices. The choice of assistive devices depends on the transfer situation and the patient's functional capability and needs.



## Functional inspection

Check the condition and function of the assistive device regularly.




## Always read the manual

Always read the manuals for all assistive devices used during a transfer.

Keep the manual where it is accessible to users of the product.

Do not leave the patient unattended during a transfer situation.


## Placement of FootStool



Place the FootStool under the patient's feet when the patient is sitting on the toilet or in a wheelchair. FootStool is recessed at the front to allow the caregiver to place a foot as close as possible to the stool. The back also has a small recess to facilitate placement close to the toilet.

## Using FootStool

FootStool complements other assistive devices such as EasyBelt/FlexiBelt, EasyGlide and SitWalk to facilitate transfers between two locations or for repositioning, for example, farther back in a chair. Since the knees are raised, the thighs are also lifted somewhat, which reduces friction and makes the transfer easier.



### Material:

Polyethylen

### Care of the product:

Clean with dishwashing liquid or disinfectant.

Wipe with a cloth.

**Note:** Do not subject the product to hard pressure or impact. This may leave indentations in the material.



Medical Device Class I

The product meets the requirements of the Medical Devices Directive 93/42/EEC.

# SystemRoMedic™

Easy-to-use assistive devices for

Easy transfers and to making life and work easier for both patients and personnel in the care sector.

Experience, innovation and training are the basis for SystemRoMedic, a total solution for every imaginable transfer situation.

**Transfer:** products for moving patients between locations.

**Positioning:** products for repositioning in the same location.

**Support:** products that provide mobility support.

**Lifting:** products adapted for lifting.

The philosophy behind SystemRoMedic is to prevent occupational injuries while improving the patient's sense of independence and dignity. Through a combination of training and a complete range of transfer-assistive devices, SystemRoMedic offers the means to improving both the work environment and the quality of care while enabling significant cost savings.

Their mission, is to help people, has always been, and will continue to be, their driving force of innovation. They love easy transfers.



[info@amilakesouthern.co.uk](mailto:info@amilakesouthern.co.uk)