

# Protective Cover

## Disposable Cover

SystemRoMedic™



**Protective Cover**



**Disposable Cover**

Max. weight for protective cover for EasySlide:  
300 kg/661 lbs

Max. weight for protective cover for RollerSlide:  
150 kg/330 lbs

Different protective covers are available for RoMedic's EasySlide and RollerSlide. Protective Cover is a multi-use cover that is made of Ultra Low Friction (ULF) material. Disposable Cover, a single-use cover that is made of plastic, can be used to advantage in e.g., emergency care or on isolation wards. The multi-use protective cover for RollerSlide is also a glide cover, since it enables the slide to glide freely.

SystemRoMedic is a series of different models of transfer assistive devices. The choice of assistive devices depends on the transfer situation and the patient's functional capability and needs.



## Functional inspection

---

### Visual inspection

Check the condition and function of the assistive device regularly. Always inspect the product after laundering. Check to ensure that seams and material are free from damage. Check to ensure that fabric is not worn or faded. If there are signs of wear, the product must be discarded.




## Always read the manual

---

Always read the manuals for all assistive devices used during a transfer.  
Keep the manual where it is accessible to users of the product.  
Do not leave the patient unattended during a transfer situation.


## Placement




Disposable Covers are on convenient rolls and are torn off one by one.



Disposable Cover and Protective Cover are pulled over EasySlide or RollerSlide (see images).



When using Disposable Cover on RollerSlide, remember that a glide cover is always needed, since it is what allows the slide to glide freely.



Note that there are different sizes of protective covers for the different products.

### Protective Cover (EasySlide)

01-1051 (01-1030,01-1020)  
01-1151 (01-1130,01-1120)  
01-1251 (01-1220,01-1230)  
01-1451 (01-1420,01-1430)  
01-1481 (01-1470,01-1480)  
01-2051 (01-2010,01-2020)  
01-2052 (01-2011,01-2021)

### Disposable Cover (EasySlide)

01-1055 (01-1030,01-1020)  
01-1155 (01-1130,01-1120)  
01-1255 (01-1220 01-1230)  
01-2055 (01-2010,01-2020)

### Sliding Cover (RollerSlide)

01-9136 (01-9120,01-9145,01-9148)  
01-9156 (01-9110,01-9155,01-9158)

### Disposable Cover (RollerSlide)

01-9138 (01-9120,01-9145 01-9148)  
01-9168 (01-9110,01-9155 01-9158)

### Care of the product:

Read the product label.



Do not use rinsing agent.

To prolong product life, avoid tumble-drying.

**NOTE!** Disposable Cover cannot be laundered, but must be discarded after use.



P O Box 4791 Ferdown Dorset BH22 9WB  
Tel: 07857 603031 Fax: 01202 877240  
Email: [info@amilakesouthern.co.uk](mailto:info@amilakesouthern.co.uk)