

Max: 150 kg/330 lbs

ReTurn7100 is used indoors for shorter transfers between wheelchair and bed, wheelchair and toilet/portable toilet/wheelchair or chair/armchair. ReTurn7100 can also be used to aid repositioning farther back in the chair. ReTurn7100 is designed for easy manoeuvring and use in confined spaces. It is easy to assemble and disassemble in just two sections for easier transport and storage.

SystemRoMedic is a series of different models of transfer assistive devices. The choice of assistive devices depends on the transfer situation and the patient's functional capability and needs.



## Functional inspection

---

### Visual inspection

Check the condition and function of the product regularly. Check to ensure that material is free from damage. If there are signs of wear, the product must be inspected and serviced.



## Always read the manual

---


Always read the manuals for all assistive devices used during a transfer.

Keep the manual where it is accessible to users of the product.

Do not leave the patient unattended during a transfer situation.

The product should be inspected yearly, and otherwise as required, by a qualified technician. Service must be performed by a qualified technician.


## Special requirements and conditions



1

- ReTurn is not intended for long transfers.
- The caregiver must always apply counter pressure when the patient is standing up, sitting down (1) and during the actual transfer (2).


**To transfer with ReTurn, the patient must be able to:**



2

- support weight while standing
- stand up and maintain balance (with support)
- grasp handles
- understand instructions


## Gripping alternatives



When several gripping alternatives are needed: The handle can be gripped horizontally at two height levels, as well as vertically between these two levels.

For tall patients: When the handle is placed at the highest position, the height of the product is raised by 20 cm. When the patient leans forward, the U-handle provides a greater area for support.


## Transport




The chassis and raising ladder can be carried separately. The chassis section can be carried with the help of the accompanying carrier strap. Pull the strap through the two holes on the underside of the base plate. Lock by pressing down on the tab.

Transport, e.g., by car, and storage are easy, since the ReTurn can be disassembled in two pieces.

## Chassis and placement of feet



The feet are placed one on either side, see foot markings on base plate.



The patient's weight must be centred on the product's support surface (all weight between the wheels); i.e., the patient must be in a standing position. The patient's weight must not be outside of the support surface.


## Brakes



There is a brake for each wheel. To brake, press down on the brake plate so that it is horizontal.

To release the brake, press on the brake plate so that it resumes its original, angled position; ensure that the wheels roll freely.

## Lower-leg support



Adjust the lower-leg support so that the upper edge of the plate is about 3 cm/two fingers below the lower edge of the kneecaps. Tighten the knob securely after adjusting.

## Transfers



**Bed to wheelchair:** Place the wheelchair at some distance from the bed to allow enough working space. Lock the wheelchair wheels. Place the patient's feet on the base plate. The patient grasps the raising ladder, leans forward, and then raises himself/herself.


Pull the ReTurn back somewhat, and then steer in towards the wheelchair.



Position the ReTurn close enough to enable good placement in the wheelchair. The caregiver applies counter pressure before the patient sits down.



**Two caregivers with ReTurnBelt:** The RoMedic ReTurnBelt allows different gripping alternatives, and helps and encourages the patient to lean forward. One caregiver helps the patient to lean forward and provides support during raising. See manual for ReTurnBelt.





**To toilet/shower chair:** For transfers to the toilet, space is allowed for the caregiver to assist the patient with the pants. The back also has a small recess that allows the patient to get close to the toilet and achieve a good seated position.




## Farther back in the chair

If the patient slides forward in the wheelchair, the ReTurn can be used to help the patient to slide farther back in the chair.



Remove/raise the wheelchair foot plates, and then run the ReTurn in towards the wheelchair. The patient places his/her feet on the foot plate, grasps the raising ladder, leans forward, and then lifts from the seat. The caregiver provides counter pressure. Slide the patient farther back in the chair when he/she has lifted slightly from the seat.

## Technical specifications



## Care of the product:

- Use a soft cloth and mild cleaning agent such as dishwashing liquid or car shampoo to clean the ReTurn. Do not use abrasive cloths or brushes to clean the base plate.
- Do not use solvents.
- To disinfect, use 70% alcohol.



P O Box 4791 Ferndown Dorset BH22 9WB

Tel: 07857 603031 Fax: 01202 877240

Email: [info@amilakesouthern.co.uk](mailto:info@amilakesouthern.co.uk)