

MiniLift125 is a mobile sit-to-stand lift which, together with a specially adapted rising vest, enables users with reduced function to rise up. The construction of MiniLift125 offers a totally unique and very comfortable sit-to-stand experience. It is also small and convenient and very easy to manoeuvre, which, in combination with the possibility to adjust the settings to the user, makes it very useful in most spaces and situations.

MiniLift125

MiniLift125 is a mobile sit-to-stand lift which has been developed to, as gently and comfortably as possible, assist the user when rising from a sitting to a standing position. The handles and the side supports offer multiple options for the user to place his/her hands. The stable anti-slip foot plate is very low and the lower leg supports as well as the lifting arm, can be adjusted to fit the user. When MiniLift125 is combined with the appropriate lifting accessories the user gets support under the feet, for the front of the lower legs and behind the back, which provides for a safe and secure sit-to-stand procedure. The construction safely moves the user forwards and upwards in a natural movement pattern and, at the same time, the leg muscles and balance is exercised. MiniLift125 is manufactured in steel for stability and strength and is suitable for users weighing up to 125 kg/275 lbs.



ACCESSORIES

MiniLift125 is used together with ThoraxSling, a specially adapted sit-to-stand vest, which is available in many different sizes. Thorax-Sling is very easy for the carer to fit and it provides support and security for the user when rising and for standing excercises. It also makes toileting chores easier. MiniLift125 and ThoraxSling provides a very gentle and pleasant sit-to-stand experience. The movement of the lift arm is designed so that the pulling motion ends comfortably in the lower back, which prevents the vest from riding up under the arms.

Disposable material for maximum hygiene and safety

ThoraxSling, with seat support is also available in a strong and durable non-woven material which is disposed of after use, when it gets soiled or when the user does not need it anymore, instead of being cleaned. This prevents the spread of multi-resistant and other types of contagious bacterial.

ThoraxSling SystemRoMedic



01-45500003 ThoraxSling, XS ThoraxSling, S 01-45500004 01-45500006 ThoraxSling, M ThoraxSling, L 01-45500007 ThoraxSling, XL 01-45500008 01-45500009 ThoraxSling, XXL SWL 300 kg / 660 lbs

Specification

Specification

ThoraxSling, with seat support

Art. No

Art. No

5 pcs

SWL 200 kg / 440 lbs



Specification 01-45600004 ThoraxSling, with seat support, S 01-45600006 ThoraxSling, with seat support, M 01-45600007 ThoraxSling, with seat support, L ThoraxSling, with seat support, XL 01-45600008 SWL 300 kg / 660 lbs

ThoraxSling, disposable





Art. No Specification 01-70200033 CalfStrap, for MiniLift125 with serial numbers from 125001653

01-45690004 ThoraxSling, with seat support, disposable, S 01-45690006 ThoraxSling, with seat support, disposable, M

01-45690007 ThoraxSling, with seat support, disposable, L

01-45690008 ThoraxSling, with seat support, disposable, XL



P O Box 4791 Ferndown Dorset BH22 9WB Tel: 07857 603031 Fax: 01202 877240 Email: info@amilakesouthern.co.uk

SystemRoMedic™ is a complete solution for easy transfers. The concept encompasses four product areas that meet most needs in different transfer situations. Transfer, between two locations. Positioning, within the same location. Support, for mobility. Lifting, both manual and mechanical.

Lenath of base

Height of footplate

Emergency lowering Weight_

Safe working load

Height _

Castor size

Lifting speed

Width of base, inner measurement

Width of base, outer measurement

CE BE BACHA AND

900 mm / 36"

980 mm / 39.2"

70 mm / 2.8"

Ø 80 mm

Art. No

445 -745 mm / 17,8" - 29,8"

590 - 890 mm / 23,6" - 35,6"

3,5 cm/s (1,4 "/s) without load

Productspecification

Manual and electrical

MiniLift125 01-400641434

_ 39 kg / 85,8 lbs 125 kg / 275 lbs