

Transfer and lifting aids

SystemRoMedic™



SystemRoMedic™ - transfer and lifting aids and accessories
June 2012



SystemRoMedic™ - simple solutions for great results.

During the 1980s, load injuries increased dramatically within the care sector. At the same time, RoMedic was founded and the very first SystemRoMedic™ transfer aid, the sliding mat EasySlide, was launched. EasySlide represented an entirely new approach and a new way of thinking within health care and when facilitating a supine transfer, and the positive market response told them they were on the right track.

25 years in the service of simplicity

dedicated to developing clever and easy-to-use assistive devices for safe and easy transfers and to making everyday life and work easier for both users and personnel in the care sector. And today, with SystemRoMedic™, it offers the market the broadest and most complete range of easy-to-use, effective and safe assistive devices for all types of patient transfers.

The challenge

Every day, care staff find themselves in situations where they are pressed for time, short of personnel and lacking knowledge and access to simple and functional assistive devices. Inadequate budgets in combination with the ambition and desire to provide a high standard of care contribute to the occurrence of various load injuries. The result is increased staff absenteeism and escalating costs. Putting an end to this negative spiral can be difficult.

The solution is simple, naturally

The transfer-assistive devices in the SystemRoMedic™ range are based on simplicity, logical thinking, ergonomic principles and design with the user and caregiver in focus. Simplicity is the key. Attention to the basic principles results in user-friendly solutions. And, in fact, it is the simple, most time-saving devices that actually get used. SystemRoMedic™ works with the laws of nature – not against them. Friction, for example, is a factor that must be considered. With smart solutions and combinations of special, functional materials, we can both reduce the friction under the user to facilitate transfer or repositioning, and increase the friction to minimize the risk of the user sliding out of bed.

Activating the user

The more strength and functional capacity the user has, the more he is able to participate in his own transfer or repositioning. And the more active the user is, the easier the transfer will be for the caregiver. The SystemRoMedic™ transfer-assistive devices help the caregiver to activate the user while stimulating a natural movement pattern. In addition to strengthening the user's muscles, this also contributes to the user's sense of independence. In this way, we can improve the quality of life for the user while at the same time making the caregiver's job easier and safer.





We love easy transfers!

The number of people requiring transfer assistance is steadily increasing. And load injuries and accidents cause huge occupational health and safety problems that account for enormous costs in the care sector.

SystemRoMedic™ offers simple solutions for great results. It's a matter of simplifying life and work and of preventing load injuries and accidents that occur in connection with patient transfers. It also has to do with activating and rehabilitating people, facilitating transfers, even for pain-sensitive users, and preventing pressure sores.

And, not least, it is a matter of improving the user's sense of independence, dignity and self-esteem.

With SystemRoMedic™, care providers can take a step in the right direction. Our products are simple, safe and renowned for their cost-effectiveness.















SystemRoMedic

SystemRoMedic™ is a complete solution that provides the simplest to the most complex scenarios, from the lightest to the heaviest. The concept encompasses products for four different categories of transfers

- **Transfer**, assistive devices for transferring users between two locations
- **Positioning**, assistive devices for repositioning of users within the same location
- **Support**, assistive devices for mobility support e.g., during sit-to-stand or gait training
- **Lifting**, assistive devices for manual and mechanical lifting of users

Find the right transfer aid for the right situation.

SystemRoMedic™ offers solutions for all transfer situations and needs. To make it easier to find and choose the right transfer aid for the right situation, please use the below matrix to combine the main products in four transfer categories with the 14 most common activities in patient transfers.

Products															
	Supine transfer	Turning in bed	Higher up in bed	Lying to sitting	Up from the floor	Support for gait training	Standing training	From bed to wheelchair	From bed to commode	From wheelchair to toilet	Rising - sitting down	Positioning in chair	In to and out of vehicles	Positioning in car	Preventing pressure sores
Transfer <i>Sliding aids for manual transfers of users between two locations.</i>															
EasyGlide	✓						✓	✓	✓				✓		
MoveMaster		✓	✓				✓	✓	✓						
EasySlide	✓	✓													✓
RollerSlide	✓														
Positioning <i>Sliding aids, turning plates and other positioning aids and accessories for manual repositioning of users within the same location.</i>															
OneWaySlide												✓	✓		
EasyGlide												✓	✓		
AutoSlide													✓	✓	
EasySlide	✓	✓	✓				✓								✓
WendyLean			✓	✓											
WendyLett		✓													✓
WendyLett2Way		✓		✓											✓
WendyLett4Way		✓	✓	✓											✓
TurnSheet		✓	✓												✓
TurnTableSoft				✓									✓		
SitWalk												✓			
ComfortCushion			✓												
AntiSlip			✓							✓					
FootStool							✓	✓	✓		✓				
TopSheet		✓													✓
		✓	✓												
Support <i>Different types of sit-to-stand aids, belts, turning plates and other support aids for mobility support e.g., during sit-to-stand or gait training.</i>															
MiniLift125, MiniLift200							✓	✓	✓	✓	✓	✓			
ThoraxSling							✓	✓	✓	✓	✓	✓			
RiseUp600EE							✓	✓	✓	✓	✓	✓			
StandUpSling							✓	✓	✓	✓	✓	✓			
QuickMove							✓	✓	✓	✓	✓	✓			
FlexiMove			✓	✓			✓	✓	✓	✓	✓	✓			
FlexiGrip				✓			✓								
StandUp					✓										
ReTurn							✓	✓	✓	✓	✓	✓			
ReTurnBelt							✓	✓	✓	✓	✓	✓			
EasyBelt/ FlexiBelt/ DisposableBelt							✓	✓	✓	✓	✓	✓			
WalkingBelt							✓	✓	✓	✓	✓	✓			



Products	Activities														
	Supine transfer	Turning in bed	Higher up in bed	Lying to sitting	Up from the floor	Support for gait training	Standing training	From bed to wheelchair	From bed to commode	From wheelchair to toilet	Rising - sitting down	Positioning in chair	In to and out of vehicles	Positioning in car	Preventing pressure sores
TurnTable								✓	✓	✓					
TurnSafe							✓	✓	✓	✓	✓				
Lifting <i>Manual lifting aids and mechanical lifting aids like ceiling lifts, rail systems, mobile lifts, lifting slings and accessories for manual or mechanical lifting of users.</i>															
LiftSeat					✓			✓	✓	✓			✓		
FlexiMove	✓	✓	✓		✓										
MiniStretcher	✓	✓	✓		✓										
RiseAtlas, stationary ceiling lift units	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
RisePorto, portable ceiling lift units	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
Castor/Pollux, freestanding lift stands	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
MilkyWay, rail system	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
MultiLift550	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
Eva, mobile lifts	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
Rino600EE	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
Vega505EE	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
HighBackSling					✓			✓		✓		✓			
ShellHBSling					✓			✓		✓		✓			
FlexibleSling					✓			✓	✓	✓		✓			
HygieneHBSling									✓	✓		✓			
BasicSling								✓				✓			
AmpSling								✓				✓			
LowBackSling								✓				✓			
HygieneSling									✓	✓					
HygieneLBSling									✓	✓		✓			
ClassicHBSling								✓				✓			
ClassicSling								✓				✓			
BariSling								✓							
StretcherSling	✓		✓		✓										
SafeHandlingSheet	✓	✓	✓	✓	✓										
TurningSling		✓													
StandingVest							✓								
WalkingVest				✓		✓	✓				✓				
BariVest						✓	✓				✓				

Sliding aids



EasyGlide 5020, 5050



EasyGlide 5060, 5070

EasyGlide

For sitting transfers from wheelchairs to beds, toilets, car seats and many other seated positions. It is especially suitable for levelling out height differences or for bridging short distances. Available in two models and several sizes.

Art. No.	Product name	Description	Unit
01-5020	EasyGlide	60 x 33 cm/ 24" x 13.2", three-part	1 Pce
01-5050	EasyGlide	75 x 33 cm/ 29.5" x 13.2", three-part	1 Pce
01-5060	EasyGlide	60 x 25 cm/ 24" x 9.8"	1 Pce
01-5070	EasyGlide	45 x 20 cm/ 17.7" x 7.8"	1 Pce

Accessories

01-4014	Storage pocket	For 5020, 5050	1 Pce
---------	----------------	----------------	-------



For supine transfers from one bed to another.

Art. No.	Product name	Description	Unit
01-5000	EasyGlide, extra long	120 x 60 cm/ 47.2" x 23.6"	1 Pce
01-5010	EasyGlide, extra long	170 x 59 cm/ 66.9" x 23.2"	1 Pce



EasyGlide

For supine transfers from one bed to another but also for positioning higher up in bed and for placing a lifting sling.

Art. No.	Product name	Description	Unit
01-5040	EasyGlide, oval	85 x 39 cm/ 33.4" x 15.3"	1 Pce
01-5042	EasyGlide, oval	85 x 39 cm/ 33.4" x 15.3"	2 Pcs



Sliding aids



MoveMaster

MoveMaster is a slightly different, open sliding mat with handles at one end. This makes it possible to help by pulling users in various transfers, such as moving them further up in bed. MoveMaster is particularly well suited for use against bare skin when transferring to shower chairs or to the toilet.

Art. No.	Product name	Description	Unit
01-2090	MoveMaster	47 x 120 cm/ 18.5" x 47.2", 3 handles	1 Pce



EasySlide

EasySlide (190 x 60) is primarily used for supine transfers, for example between beds and x-ray or operating tables, shower trolleys, stretchers and similar. EasySlide is also a very effective aid for turning users in bed or on x-ray or operating tables. Disposable covers for EasySlide provide for improved hygiene and less laundry.



EasySlide Disposable



Protective Cover



Disposable cover

Art. No.	Product name	Description	Unit
01-1020	EasySlide	190 x 60 cm/ 76" x 24"	1 Pce
01-1030	EasySlide (1020 + 1051)	190 x 60 cm/ 76" x 24" complete	1 Set
01-2040	EasySlide Disposable	50 x 70 cm	2 x 60 pcs
01-2045	EasySlide Disposable	50 x 70 cm inc. wall/table dispenser	2 x 60 pcs

Accessories

01-1051	Protective cover	For 1020, 1030	1 Pce
01-1055	Disposable cover	For 1020, 1030	2 x 50 Pcs
01-1151	Protective cover	For 1120*, 1130*	1 Pce
01-1155	Disposable cover	For 1120*, 1130*	2 x 90 Pcs
01-1085	WallRack	For storing of sliding mats and disposable covers for EasySlide and RollerSlide	1 Pce



* 1120, 1130, EasySlide 100 x 60, discontinued



WallRack


Sliding aids



RollerSlide 01-9145



RollerSlide, flexi



Disposable cover



WallRack

RollerSlide

For supine transfers of inert users between a bed and an x-ray or operating table, shower trolley, ambulance stretcher or similar. Bridges gaps and moves with the users during the transfer.

RollerSlide consists of a padded sliding board with a cover made from ULF (Ultra Low Friction), a sophisticated sliding material that slides round the board during transfer. Available in several versions and sizes. RollerSlide, flexi is a lighter and bendable version of RollerSlide.

RollerSlide, light is a light-weight version of RollerSlide, flexi.

Art. No.	Product name	Description	Unit
01-9145	RollerSlide	175 x 50 cm/70 x 20", foldable	1 Pce
01-9155	RollerSlide	85 x 50 cm/ 34" x 20"	1 Pce
01-9148	RollerSlide, flexi	175 x 50 cm/ 70" x 20", foldable	1 Pce
01-9158	RollerSlide, flexi	85 x 50 cm/ 34" x 20"	1 Pce
01-9120	RollerSlide, light	175 x 50 cm/ 70" x 20", foldable	1 Pce
01-9110	RollerSlide, light	85 x 50 cm/ 34" x 20"	1 Pce

Accessories

01-9136	Sliding cover	For 9120, 9145, 9148	1 Pce
01-9156	Sliding cover	For 9110, 9155, 9158	1 Pce
01-9138	Disposable cover	For 9120, 9145, 9148	2 x 50 Pcs
01-9168	Disposable cover	For 9110, 9155, 9158	2 x 115 Pcs
01-1085	WallRack	For storing of sliding mats and disposable covers for EasySlide and RollerSlide	1 Pce





Sliding aids



OneWaySlide 01-2075



OneWaySlide 01-2069

OneWaySlide

OneWaySlide has many functions. The main one improves sitting positions in chairs or wheel-chairs. The sliding mat effectively retards forward movement yet slides backwards, making it easy for the carer to help the user move into a more comfortable position in the chair. Also available in a tube-shaped design without handles for the more active user.

Art. No.	Product name	Description	Unit
01-2074	OneWaySlide	110 x 30 cm / 44" x 12", 6 handles	1 Pce
01-2075	OneWaySlide	130 x 40 cm / 52" x 16", 6 handles	1 Pce
01-2076	OneWaySlide	130 x 45 cm / 52" x 18", 6 handles	1 Pce
01-2077	OneWaySlide	130 x 55 cm / 52" x 22", 6 handles	1 Pce
01-2069	OneWaySlide, tubular	50 x 40 cm / 20" x 16"	1 Pc
01-2070	OneWaySlide, tubular	50 x 45 cm / 20" x 18"	1 Pce



EasyGlide

EasyGlide are used in pairs to facilitate positioning further back in chair. The smaller size makes it easier to place the sliding boards under the user's thighs and seat.

Art. No.	Product name	Description	Unit
01-5045	EasyGlide, oval	45 x 20 cm / 17.7" x 7.9"	2 Pcs



AutoSlide

Facilitates getting in and out of cars, allowing its users to easier turn and move sideways. The seat cushions are firmly attached to the backs of car seats, so users cannot slide forward.

Art. No.	Product name	Description	Unit
01-8100	AutoSlide	48 x 48 cm / 19.2" x 19.2"	1 Pce



Sliding aids



EasySlide

EasySlide (100 x 70) is mainly designed for users requiring only some care who need to turn and move in bed. It performs best on a rough surface. EasySlide reduces the friction on both shoulder and hip, which in turn reduces the risk of pressure sores. Thanks to these features, EasySlide is very gentle for people in pain but also for the carer, thus avoiding the risks associated with lifting.



Disposable cover



WallRack

Art. No.	Product name	Description	Unit
01-1220	EasySlide	100 x 70 cm/ 40" x 28"	1 Pce
01-1230	EasySlide (01-1220 + 01-1251)	100 x 70 cm/ 40" x 28", complete	1 Set

Accessories

01-1251	Protective cover	For 1220, 1230	1 Pce
01-1255	Disposable cover	For 1220, 1230	2 x 90 Pcs
01-1085	WallRack	For storing of sliding mats and disposable covers for EasySlide and RollerSlide	1 Pce



EasySlide 01-1430

EasySlide

EasySlide (140 x 90) is designed for normal-weight users with limited to severely reduced abilities to participate in the turning process. These people are often permanently confined to the bed, and it is therefore important to prevent pressure sores. EasySlide (200 x 90) has extra padding in order to further reduce the risks of pressure sores. The padding distributes the pressure over the entire body and reduces pressure, friction and the risk of shearing on the especially vulnerable pressure points: the shoulders, hips, head, and feet.



EasySlide 01-1470

Art. No.	Product name	Description	Unit
01-1420	EasySlide	140 x 90 cm/ 56" x 36", 2 fastening straps	1 Pce
01-1430	EasySlide (1420 + 1451)	140 x 90 cm/ 56" x 36", complete	1 Set
01-1470	EasySlide	200 x 90 cm/ 80" x 36", 4 fastening straps	1 Pce
01-1480	EasySlide (1470 + 1481)	200 x 90 cm/ 80" x 36", complete	1 Set

Accessories

01-1451	Protective cover	For 1420, 1430	1 Pce
01-1481	Protective cover	For 1470, 1480	1 Pce
01-1599	LockingStrap	Fastening straps, sparepart for 1420, 1430, 1470, 1480	1 Pair
01-1085	WallRack	For storing of sliding mats and disposable covers for EasySlide and RollerSlide	1 Pce



WallRack

Sliding aids



EasySlide 01- 2020



EasySlide 01- 2015



EasySlide disposable mats 01-2040



Disposable cover



WallRack

EasySlide

EasySlide (50 x 60), (50 x 70) and (25 x 60) are excellent aids for turning in bed, for moving users higher up in bed or for transferring them from the bed to the wheelchair. EasySlide is a small, practical aid that encourages activation. For the carer, EasySlide decreases the number of assisted positionings and eliminates unnecessary lifting, which spares both back and shoulders.

Art. No.	Product name	Description	Unit
01-2010	EasySlide	50 x 60 cm/ 19.7" x 23.6", polyester	1 Pce
01-2020	EasySlide	50 x 60 cm/ 19.7" x 23.6", microfiber	1 Pce
01-2030	EasySlide (2020 + 2051)	50 x 60 cm/ 19.7" x 23.6", complete	1 Set
01-2011	EasySlide	50 x 70 cm/ 19.7" x 27.6", polyester	1 Pce
01-2021	EasySlide	50 x 70 cm/ 19.7" x 27.6", microfiber	1 Pce
01-2015	EasySlide	25 x 60 cm/ 9.8" x 23.6", polyester	1 Pce
01-2040	EasySlide (disposable)	50 x 70cm/19.7 x 27.6", Polyeten PE	2 x 60pcs
01-2045	EasySlide (disposable) with dispenser	50 x 70cm/19.7 x 27.6" Polyeten PE in steel for wall or tabletop.	2 x 60 pcs
Accessories			
01-2051	Protective cover	For 01-2010, 01-2020, 01-2030	1 Pce
01-2052	Protective cover	For 01- 2011, 01-2021	1 Pce
01-2055	Disposable cover	For 01-2010, 01- 2020, 01- 2030	100 Pcs
01-1085	WallRack	For storing of sliding mats and disposable covers for EasySlide and RollerSlide	1 Pce



WendyLean

The patented, self-retracting handles are very handy when you need to help a user to sit up in bed. The underside of the pillowcase is made of smooth satin which makes it very easy to pull the pillow in under the shoulders for raising the upper body of the user.



Art. No.	Product name	Description	Unit
01-1680	WendyLean	50 x 60 cm/ 20" x 24", grey/white striped	1 Pce




Sliding aids

WendyLett

WendyLett base sheet is a satin sliding sheet. The sliding takes place sideways along the direction of the stripes. WendyLett is intended for users who have difficulties turning in bed.

WendyLett has a smooth, silky central part, enabling turning with minimal effort. Sliding sideways is made easy without risk of sliding off the bed by the anti-slip material at both sides of the sheet. The material also prevents the user from sliding down in bed. WendyLett is bedded in the bed like a normal sheet with the smooth surface in the middle of the bed.

For users who need more assistance, WendyLett can be combined with WendyLett2Way draw-sheet or WendyLett4Way sliding drawsheet.



WendyLett 01-1641, 01-1651



Art. No.	Product name	Description	Unit
01-1609	WendyLett	Sliding sheet, 100 x 200 cm/ 40" x 80", white, sliding surface 60 cm/ 24", fits 80-90 cm/ 32"-36" bed	1 Pce
01-1614	WendyLett	Sliding sheet, 140 x 200 cm/ 56" x 80", white, sliding surface 60 cm/ 24", fits 80-90 cm/ 32"-36" bed	1 Pce
01-1629	WendyLett	Sliding sheet, 100 x 200 cm/ 40" x 80", grey/white striped, sliding surface 60 cm/ 24", fits 80-90 cm/ 32"-36" bed	1 Pce
01-1634	WendyLett	Sliding sheet, 140 x 200 cm/ 56" x 80", grey/white striped, sliding surface 60 cm/ 24", fits 80-90 cm/ 32"-36" bed	1 Pce
01-1639	WendyLett	Sliding sheet, 200 x 200 cm/ 80" x 80", grey/white striped, sliding surface 70 cm/ 28", fits 105 cm/ 42" bed	1 Pce
01-1667	WendyLett	Sliding sheet, 100 x 220 cm/ 40" x 87", grey/white striped, sliding surface 100 cm/ 40", fits 120 cm/ 48" bed	1 Pce
01-1668	WendyLett	Sliding sheet, 140 x 220 cm/ 56" x 87", grey/white striped, sliding surface 100 cm/ 40", fits 120 cm/ 48" bed	1 Pce
01-1669	WendyLett	Sliding sheet, 200 x 220 cm/ 80" x 87", grey/white striped, sliding surface 100 cm/ 40", fits 120 cm/ 48" bed	1 Pce
01-1641	WendyLett	Sliding sheet, 90 x 200 cm/ 36" x 80", grey/white striped, tailored w. drawstring, sliding surface 70 cm/ 28", fits 90 cm/ 36" bed	1 Pce
01-1651	WendyLett	Sliding sheet, 105 x 200/ 42" x 80", grey/white striped, tailored w. drawstring, sliding surface 80 cm/ 32", fits 105 cm/ 42" bed	1 Pce



PREVENTS
pressure sores

Sliding aids




WendyLett2Way

WendyLett2Way is a drawsheet made of a material which is adapted to WendyLett base sheet so that sliding down in bed is much reduced. Sliding takes place sideways, along the direction of the stripes. WendyLett2Way is placed on top of a WendyLett base sheet. The combination makes it easy for the caregiver to turn the user sideways.

Art. No.	Product name	Description	Unit
01-1620	WendyLett2Way	Drawsheet 90 x 180 cm/ 36" x 72", white	1 Pce
01-1631	WendyLett2Way	Drawsheet 100 x 200 cm/ 40" x 80", white/grey striped	1 Pce
01-1636	WendyLett2Way	Drawsheet 140 x 200 cm/ 56" x 80", white/grey striped	1 Pce

Accessories

01-70200007	RoMedicClips	Accessories to WendyLett and TurnSheet	1 Pair
-------------	--------------	--	--------



WendyLett4Way

WendyLett4Way is a sliding drawsheet to be combined with WendyLett base sheet. WendyLett4Way offers multi-directional sliding, as indicated by the checked pattern. WendyLett4Way is placed on top of a WendyLett base sheet. The combination allows for multi-directional moves and facilitates procedures like turning in bed, moving higher up/lower in bed or other re-adjusting of a lying user, even with heavier or passive users.

Art. No.	Product name	Description	Unit
01-1637	WendyLett4Way	Sliding drawsheet, 140 x 200 cm/ 56" x 80", grey/white checked	1 Pce
01-1647	WendyLett4Way	Sliding drawsheet, 200 x 200 cm/ 80" x 80", grey/white checked	1 Pce
01-1648	WendyLett4Way, inco	Sliding drawsheet/incontinence protection with extra fabric on both sides, 75 x 200 cm/ 29.5" x 80", grey/white checked	1 Pce
01-1649	WendyLett4Way, inco with handles	Sliding drawsheet/incontinence protection with 7 handles on each side, 140 x 120 cm/ 56" x 47.2", grey	1 Pce
01-1657	WendyLett4Way, with handles	Sliding drawsheet with 4 handles on each side, 140 x 190 cm/ 56" x 75", grey/white checked	1 Pce

Accessories

01-70200007	RoMedicClips	Accessories to WendyLett and TurnSheet	1 Pair
-------------	--------------	--	--------



Sliding aids



TurnSheet

TurnSheet (65 x 130) is used to enable turning and changing position in bed. This is a very popular product, especially for people with back problems and for pregnant women with pelvic conditions.

Art. No.	Product name	Description	Unit
01-2060	TurnSheet	65 x 130 cm/ 26" x 52"	1 Pce
01-2080	TurnSheet (01-2060 + 01- 2081)	65 x 130 cm/ 26" x 52", complete	1 Set

Accessories

01-2081	Mattress cover	For 2060, 2080, L 60 cm/ 23.1", fits 90 cm/ 35.4" bed	1 Pce
---------	----------------	---	-------



TurnSheet

TurnSheet (150 x 200) is developed for heavier users, people with severely reduced functional abilities and people at risks for pressure sores who need assistance with turning in bed.

Art. No.	Product name	Description	Unit
01-1592	TurnSheet	150 x 200 cm/ 60" x 80"	1 Pce
01-1580	TurnSheet (01-1592 + 01-1583)	150 x 200 cm/ 60" x 80", complete	1 Set

Accessories

01-1583	Mattress cover	For 1592, 1580, L 140 cm/ 56", fits 90 cm/ 35.4" bed	1 Pce
01-70200007	RoMedicClips	Accessories to WendyLett and TurnSheet	1 Pair



RoMedicClips are used as an aid for turning users. RoMedicClips can be attached to most sling bars and can be used on SystemRoMedic's drawsheets. The carer can perform the turn without any load.



Sliding aids



TurnSheet

A pain-sensitive user, requiring intensive and frequent care can be a challenge for carers. With TurnSheet (200 x 150) little strength is needed and unnecessary lifting is avoided in turning. This device is often used in hospitals and in home care.

Art. No.	Product name	Description	Unit
01-1592	TurnSheet	200 x 150 cm/ 80" x 60"	1 Pce
01-1590	TurnSheet (01-1592 + 01-1593)	200 x 150 cm/ 80" x 60", complete	1 Set
01-1589	TurnSheet	200 x 240 cm/ 80" x 94.4"	1 Pce

Accessories

01-1593	Mattress cover	For 1592, 1590 and 1510, 1515 (TopSheet), L 200 cm/ 80", fits 90 cm / 35.4" bed	1 Pce
01-1594	Mattress cover	For 1592, 1589, and 1510, 1515 (TopSheet), L 200 cm/ 80", fits 105 cm / 41.3" bed	1 Pce
01-1595	Mattress cover	For 1592 and 1589, L 200 cm/ 80", fits 120 cm / 47.2" bed	1 Pce
01-70200007	RoMedicClips	Accessories to WendyLett and TurnSheet	1 Pair



RoMedicClips are used as an aid for turning users. RoMedicClips can be attached to most sling bars and can be used on SystemRoMedic's drawsheets. The carer can perform the turn without any load.



Turning plates



TurnTableSoft

TurnTableSoft makes it easier to turn when sitting down. TurnTableSoft is constructed of two round discs with a friction-reducing material on the inside, enabling the discs to rotate smoothly over each other.

Art. No.	Product name	Description	Unit
01-5086	TurnTableSoft	Ø 45 cm/ Ø 18"	1 Pce
01-5087	TurnTableSoft	Ø 40 cm/ Ø 16"	1 Pce



Other positioning aids



SitWalk 01-6080

SitWalk

This device is named after the method of moving a user backwards in the chair to a more comfortable sitting position. Lean the user sideways and move one buttock at a time, facilitating a hip-hitch.

Art. No.	Product name	Description	Unit
01-6080	SitWalk	140 x 13 cm/ 56" x 5.2", 8 handles	1 Pce
01-6085	SitWalk, extra long	245 x 12 cm/ 96.4" x 4.7", 8 handles, accessory for StandUp	1 Pce



ComfortCushion

With ComfortCushion under the thighs, the user maintains a relaxed posture even when the head end of the bed is raised. The cushion also prevents the user from sliding down in bed.

Art. No.	Product name	Description	Unit
01-8050	ComfortCushion	40 x 50 cm/ 16" x 20"	1 Pce



Other positioning aids



AntiSlip

AntiSlip comes in handy when you need a good grip for moving higher up in bed. It can also be used as a restraining surface between materials in many of our sliding products.

Art. No.	Product name	Description	Unit
01-8040	AntiSlip	45 x 38 cm/ 18" x 15.2"	1 Pce
01-8045	AntiSlip	45 x 38 cm/ 18" x 15.2"	10 Pcs
01-8046	AntiSlip	45 x 38 cm/ 18" x 15.2"	30 Pcs



FootStool 01-8090



FootStool 01-8095

FootStool

FootStool is used to raise the position of the user's knees to enable an easier and safer transfer in a seated position. FootStool also provides excellent support under the user's feet during visits to the toilet or when repositioning further back in a chair or wheelchair.

Art. No.	Product name	Description	Unit
01-8090	FootStool	32 x 38 x 10 cm/ 12.8" x 15.2" x 4"	1 Pce
01-8095	FootStool	32 x 38 x 5 cm/ 12.8" x 15.2" x 2"	1 Pce



Other positioning aids



TopSheet

TopSheet is used for incontinent, bedridden users who must be turned often. TopSheet is also excellent when you need to hold a heavy user in a lateral position, for instance when assisting with personal hygiene or dressing wounds.


Art. No.	Product name	Description	Unit
01-1510	TopSheet	200 x 160 cm/ 80" x 64", ULF on the underside	1 Pce
01-1515	TopSheet	155 x 140 cm/ 61" x 56", ULF on the underside	1 Pce
01-1520	TopSheet	144 x 140 cm/ 57.6" x 56"	1 Pce
01-1511	TopSheet (01-1510 + 01-1593)	200 x 160 cm/ 80" x 64", complete with Mattress cover (TurnSheet)	1 Set

Accessories

01-1509	Lower part for TopSheet 1510	140 x 160 cm/ 56" x 64"	1 Pce
01-1514	Lower part for TopSheet 1515	85 x 140 cm/ 33" x 56"	1 Pce
01-1519	Lower part for TopSheet 1520	80 x 140 cm/ 31" x 56"	1 Pce



Accessories



IncoSheet

IncoSheet is a washable incontinence protection to be placed on the bed and over existing bed equipment. IncoSheet is available in several sizes, including a size that fits a wheelchair cushion. The material is highly absorbent and distributes fluids over a large area.

Art. No.	Product name	Description	Unit
01-1530	IncoSheet	75 x 90 cm/ 30" x 36"	1 Pce
01-1531	IncoSheet	75 x 200 cm/ 30" x 80", extra fabric on both sides	1 Pce
01-1532	IncoSheet	75 x 150 cm/ 30" x 60"	1 Pce
01-1533	IncoSheet	40 x 83 cm/ 16" x 33.2"	1 Pce



Sit-to-stand aids



SlingBarWrap MiniLift



CalfStrap

MiniLift200

MiniLift200 is a mobile sit-to-stand lift developed for helping users, as gently as possible, to rise from a seated to a standing position. Used in combination with the right accessories, the user gets support under the feet, for knees and lower legs and behind the back. The bar of MiniLift200 can be provided with a soft protection, SlingBarWrap MiniLift, which is sometimes needed for users with uncontrolled movements. SlingBarWrap MiniLift is fastened with velcro and can easily be removed for cleaning.

Art. No.	Product name	Description	Unit
01-401100334	MiniLift200	Max 200 kg/ 440 lbs	1 Pce

Accessories

See page 21	ThoraxSling	Sit-to-stand vest	
01-70200012	SlingBarWrap MiniLift	For MiniLift200	1 Pce
01-70200033	CalfStrap	For MiniLift200 as from serial number 200001277	1 Pce
01-70200029	External charger Linak	For MiniLift200/ MiniLift125	1 Pce



CalfStrap

MiniLift125

MiniLift125 is a mobile sit-to-stand lift developed for helping users, as gently as possible, to rise from a seated to a standing position. Used in combination with the right accessories, the user gets support under the feet, for knees and lower legs and behind the back.

Art. No.	Product name	Description	Unit
01-400641434	MiniLift125	Max 125 kg/ 275 lbs	1 Pce

Accessories

See page	ThoraxSling	Sit-to-stand vest	
01-70200033	CalfStrap	For MiniLift125 as from serial number 125001653	1 Pce
01-70200029	External charger Linak	For MiniLift125/ MiniLift200	1 Pce



Sit-to-stand aids



ThoraxSling

ThoraxSling

ThoraxSling is a sit-to-stand vest for MiniLift125 and MiniLift200 that is easy for the caregiver to apply, while giving the user an extra feeling of security. It is designed to meet the need for support and security during rising and to facilitate visits to the toilet. ThoraxSling is available in two models; one that provides support behind the back and one that also gives support under the seat, and several sizes, all for users weighing max 300 kg/660 lbs.

ThoraxSling, with seat support is also available in a disposable material which is disposed after use instead of being cleaned. This prevents the spread of multi-resistant and other types of contagious bacteria. ThoraxSling with seat support disposable is intended for short term use and for users weighing max 200 kg/440 lbs.

Art. No.	Product name	Description	Unit
01-45500003	ThoraxSling	Polyester, X-Small	1 Pce
01-45500004	ThoraxSling	Polyester, Small	1 Pce
01-45500006	ThoraxSling	Polyester, Medium	1 Pce
01-45500007	ThoraxSling	Polyester, Large	1 Pce
01-45500008	ThoraxSling	Polyester, X-Large	1 Pce
01-45500009	ThoraxSling	Polyester, XX-Large	1 Pce
01-45600004	ThoraxSling	With seat support, polyester, Small	1 Pce
01-45600006	ThoraxSling	With seat support, polyester, Medium	1 Pce
01-45600007	ThoraxSling	With seat support, polyester, Large	1 Pce
01-45600008	ThoraxSling	With seat support, polyester, X-Large	1 Pce
01-45690004	ThoraxSling	With seat support, disposable, Small	5 Pcs
01-45690006	ThoraxSling	With seat support, disposable, Medium	5 Pcs
01-45690007	ThoraxSling	With seat support, disposable, Large	5 Pcs
01-45690008	ThoraxSling	With seat support, disposable, X-Large	5 Pcs



ThoraxSling,
with seat support



ThoraxSling,
with seat support,
disposable



Sit-to-stand aids



RiseUp600EE

RiseUp600EE is a mobile sit-to-stand lift developed for helping the heavier user, as gently and comfortably as possible, to rise from a seated to a standing position. Together with the sit-to-stand vest StandUpSling even a user with much reduced ability can have the satisfaction of being able to stand.

Art. No.	Product name	Description	Unit
01-60500001	RiseUp600EE	Max 270 kg/ 600 lbs	1 Pce

Accessories

01-45200004 -- 9	StandUpSling	Sit-to-stand vest	1 Pce
------------------	--------------	-------------------	-------



StandUpSling

StandUpSling is a sit-to-stand vest that is easy for the caregiver to apply, while giving the user an extra feeling of security. It is designed to provide support behind the back during rising with a sit-to-stand type lift, e.g. RiseUp600EE, and to facilitate visits to the toilet. StandUpSling enables a gentle rising in spite of reduced function. However, the user must have some core/trunk stability. StandUpSling is available in several sizes, all for users weighing max 300kg/660 lbs.



StandUpSling Padding can be applied to the StandUpSling. It gives StandUpSling an extra layer in soft material that attaches to the vest and is easy for the caregiver to apply. StandUpSling Padding is used mainly for users who are sensitive to pressure.

Art. No.	Product name	Description	Unit
01-45200004	StandUpSling	Polyester, Small	1 Pce
01-45200006	StandUpSling	Polyester, Medium	1 Pce
01-45200007	StandUpSling	Polyester, Large	1 Pce
01-45200008	StandUpSling	Polyester, X-Large	1 Pce
01-45200009	StandUpSling	Polyester, XX-Large	1 Pce

Accessories

01-45273306	StandUpSling Padding	Small/ Medium	1 Pce
01-45273307	StandUpSling Padding	Large/ X-Large	1 Pce



Sit-to-stand aids



QuickMove

QuickMove is a safe and easy-to-use aid for sit-to-stand and transfer of users with impaired balance and reduced lower limb strength. Based upon natural movement patterns, QuickMove supports the user when rising, using his/her own strength, and then transfers the user, either standing up or sitting down, from bed to chair or wheelchair or to and from the toilet. The user must have some strength in the arms and be able to weightbear to actively participate in the procedure.

BackSupport is an accessory which can be used as torso/side support for QuickMove.

BackSupport provides extra stability and security for users during transfers with QuickMove.

Another alternative is SupportStraps which, together with EasyBelt/FlexiBelt/FlexiMove, can be used as torso/side support.

Art. No.	Product name	Description	Unit
01-400801334	QuickMove	Max 170 kg/ 375 lbs	1 Pce

Accessories

01-70200037	BackSupport	Torso/side support for QuickMove	1 Pce
01-6012 -- 01-6015	EasyBelt	Extra support for QuickMove	1 Pce
01-6023 --01-6026	FlexiBelt	Extra support for QuickMove	1 Pce
01-6072 -- 01-6076	FlexiBelt, hug	Extra support for QuickMove	1 Pce
01-3012, 01-3020,01-3030	Fleximove	Extra support for QuickMove	1 Pce
01-6061	SupportStraps	Accessory for belts	1 Pair



BackSupport



FlexiMove

FlexiMove has multiple areas of application. For example as support when helping someone get up from a bed or chair, when moving higher up in bed or when transferring to and from a wheelchair – in that case, often on a sliding board. FlexiMove is available in two versions with different inside material.

Art. No.	Product name	Description	Unit
01-3012	FlexiMove	20 x 60 cm/ 8" x 24", nylon	1 Pce
01-3020	FlexiMove	20 x 50 cm/ 8" x 20", polyester velvet	1 Pce
01-3030	FlexiMove	20 x 60 cm/ 8" x 24", polyester velvet	1 Pce

Accessories

01-6061	SupportStraps	Accessory for belts and FlexiMove together with QuickMove	1 Pair
---------	---------------	---	--------



Sit-to-stand aids



FlexiGrip

FlexiGrip is the solution for users who have physical difficulties getting out of bed. Step by step they gradually pull themselves up from a lying to a sitting position. FlexiGrip can also be used to assist in transfers from wheelchair to bed.

Art. No.	Product name	Description	Unit
01-6090	FlexiGrip	175 x 9 cm/ 70" x 3.6", 8 handles	1 Pce

Accessories

01-6091	Extension strap	For FlexiGrip	1 Pce
---------	-----------------	---------------	-------



StandUp

StandUp is a lift-chair and a brand new kind of assistive device in the SystemRoMedic range. Using StandUp, a single caregiver can quickly, safely and easily assist a person lying on the floor, unable to get up on their own. In a secure and comfortable way, StandUp lifts the user up into a seated position, ready for rising or transferring to a wheelchair. StandUp is both compact and light-weight. It weighs only 11,5 kg/25 lbs, which is less than half compared to a portable mobile lift. At the same time, it is strong and sturdy enough to manage lifting users weighing up to 160 kg/352 lbs. StandUp is delivered complete with carrying strap and an extra long version of SitWalk (01- 6085) for helping the user get up on the seat.

Art. No.	Product name	Description	Unit
01-60800001	StandUp	45 x 74 x 19 cm/ 17,7" x 29,1" x 7,4", folded	1 Pce

Accessories

01-0951	SpaciousBag	Large bag with wheels for convenient transportation and storage of StandUp and other assistive devices	1 Pce
---------	-------------	--	-------



Sit-to-stand aids



ReTurn 01-7100 with U-handle



ReTurn 01-7500



ReTurn 01-7400



HeelStrap

ReTurn

ReTurn is a multi-functional, pivoting and transportation platform. It is also a rehabilitation device. For users who still retain some strength, ReTurn stimulates the natural movement pattern when progressing from sit-to-stand and, at the same time, it also strengthens a range of muscles. ReTurn provides a secure and active transfer when the user has to stand up and be moved to a bed, wheelchair or toilet. It can also be used to aid repositioning further back in the chair. For the carer, it is an easy procedure which eliminates all heavy lifting procedures. ReTurn is available in three models, all of them designed for easy manoeuvring and use in confined spaces. They are also easy to assemble and disassemble for easier transport and storage.

- ReTurn 01-7100 is excellent for the taller user or when the user stands up from a high bed or chair.
- ReTurn 01-7500 has a standard height rising ladder that suits most adults.
- ReTurn 01-7400 has a lower rising ladder which suits the shorter adults and children. It also works well in special care situations, such as in x-ray departments.

ReTurn 01-7500 and ReTurn 01-7400 have the same base plate which makes it possible to change between low and standard height rising ladder for maximum utilization.

All models of ReTurn are for users weighing max 150 kg/ 330 lbs.

Art. No.	Product name	Description	Unit
01-7100	ReTurn 01-7100	Transfer platform, grey, complete	1 Pce
01-7500	ReTurn 01-7500	Transfer platform, red, complete with standard height rising ladder	1 Pce
01-7400	ReTurn 01-7400	Transfer platform, red, complete with low rising ladder	1 Pce
01-7501	Rising ladder	Standard height; fits ReTurn7500 from serial no. 373317697500500001 and ReTurn7400 from serial no. 373317697400500001	1 Pce
01-7401	Rising ladder	Low; fits ReTurn7500 from serial no.373317697500500001 and ReTurn7400 from serial no.373317697400500001	1 Pce

Accessories

01-6031 - 01- 6036	ReTurnBelt	For 01-7100, 01-7400, 01-7500	1 Pce
01-6061	SupportStraps	Accessory for belts	1 Pair
01-7215	Padding for leg support	For 01-7100, 01-7400, 01-7500 and TurnSafe	1 Pair
01-7250	HeelStrap	For 01-7400, 01-7500	1 Pce
01-0954	Bag	For 01-7400, 01-7500	1 Pce
01-7225	U-handle	For 01-7100	1 Pce



Belts



ReTurnBelt

ReTurnBelt is used for sit-to-stand, transfers and gait training procedures. ReTurnBelt is an excellent complement to ReTurn when more support and stability is needed. ReTurnBelt has an anti-slip inside of PU-coated nylon and can therefore be used in environments where moisture and dampness are present or when there are special hygiene requirements.

Art. No.	Product name	Description	Unit
01-6031	ReTurnBelt	XXS, 88 cm/ 34.6", 2 x 2 handles	1 Pce
01-6032	ReTurnBelt	XS, 100 cm/ 39.4", 2 x 2 handles	1 Pce
01-6033	ReTurnBelt	S, 115 cm/ 45.3", 2 x 2 handles	1 Pce
01-6034	ReTurnBelt	M, 150 cm/ 59.1", 2 x 2 handles	1 Pce
01-6035	ReTurnBelt	L, 155 cm/ 61", 2 x 2 handles	1 Pce
01-6036	ReTurnBelt	XL, 180 cm/ 70.9", 2 x 2 handles	1 Pce

Accessories

01-6061	SupportStraps	Accessory for belts	1 Pair
---------	---------------	---------------------	--------



EasyBelt

EasyBelt is used for sit-to-stand transfers or for moving from the wheelchair to the bed, toilet or car. EasyBelt has an anti-slip inside of PU-coated nylon and can therefore be used in environments where moisture and dampness are present or when there are special hygiene requirements.

Art. No.	Product name	Description	Unit
01-6012	EasyBelt	XS, 75 cm/ 29.5", 5 handles	1 Pce
01-6013	EasyBelt	S, 105 cm/ 41.3", 7 handles	1 Pce
01-6014	EasyBelt	M, 120 cm/ 47.2", 9 handles	1 Pce
01-6015	EasyBelt	L, 135 cm/ 53.1", 11 handles	1 Pce
01-6002	EasyBelt, hug	XS, 75 cm/ 29.5", 3 + 2 movable handles	1 st
01-6003	EasyBelt, hug	S, 105 cm/ 41.3", 4 + 3 movable handles	1 st
01-6004	EasyBelt, hug	M, 120 cm/ 47.2", 5 + 4 movable handles	1 st
01-6005	EasyBelt, hug	L, 135 cm/ 53.1", 6 + 5 movable handles	1 st
01-6006	EasyBelt, hug	XL, 150 cm/ 59.1", 7 + 6 movable handles	1 st

Accessories

01-6061	SupportStraps	Accessory for belts	1 Pair
---------	---------------	---------------------	--------



Belts



FlexiBelt

FlexiBelt is used for sit-to-stand transfers or for moving from the wheelchair to the bed, toilet or car. FlexiBelt has a soft innerside of polyester velvet which makes it very comfortable to wear.

FlexiBelt, hug is a new version of SystemRoMedic's classic FlexiBelt. FlexiBelt, hug features a new, patented design with horizontal handles that slide freely in loops to make the belt "hug" on to the user instead of sliding up during transfer situations. With FlexiBelt, hug, the carer can provide even more support without using more effort. The new design actually does all the work by itself.

Art. No.	Product name	Description	Unit
01-6023	FlexiBelt	S, 105 cm/ 41.3", 7 handles	1 Pce
01-6024	FlexiBelt	M, 120 cm/ 47.2", 9 handles	1 Pce
01-6025	FlexiBelt	L, 135 cm/ 53.1", 11 handles	1 Pce
01-6026	FlexiBelt	XL, 150 cm/ 59.1", 13 handles	1 Pce
01-6072	FlexiBelt, hug	XS, 75 cm/ 29.5", 3 + 2 movable handles	1 Pce
01-6073	FlexiBelt, hug	S, 105 cm/ 41.3", 4 + 3 movable handles	1 Pce
01-6074	FlexiBelt, hug	M, 120 cm/ 47.2", 5 + 4 movable handles	1 Pce
01-6075	FlexiBelt, hug	L, 135 cm/ 53.1", 6 + 5 movable handles	1 Pce
01-6076	FlexiBelt, hug	XL, 150 cm/ 59.1", 7 + 6 movable handles	1 Pce



FlexiBelt, hug

Accessories

01-6061	SupportStraps	Accessory for belts	1 Pair
---------	---------------	---------------------	--------



SupportStraps

SupportStraps is a multi-functional accessory providing extra support and security for users during standing, gait training and transfers with SystemRoMedic's belts, sit-to-stand aids and transfer platforms. SupportStraps can be used as:

- leg harness for EasyBelt/FlexiBelt,
- torso/side support when using QuickMove together with EasyBelt/FlexiBelt/FlexiMove
- side support when using ReTurn together with ReTurnBelt



Art. No.	Product name	Description	Unit
01-6061	SupportStraps	Accessory for belts	1 Pair



Belts



WalkingBelt

WalkingBelt is used to support children with impaired strength and stability during gait training. WalkingBelt, which provides support around the torso, is held in place with a padded leg harness. The length of the two padded handles can be adapted to give the caregiver a comfortable working position when walking close to the child holding the handles to provide the child with additional support and stability during gait training. WalkingBelt is not to be used as an assistive device for lifting. WalkingBelt has a soft polyester velvet lining that feels pleasant against the skin and makes it very comfortable to wear. The material is also anti-slip treated to stay in position and prevent slipping.

Art. No.	Product name	Description	Unit
01-6041	WalkingBelt	XXS, 60 cm/ 23,6"	1 Pce
01-6042	WalkingBelt	XS, 70 cm/ 27,6"	1 Pce



DisposableBelt

DisposableBelt is intended for short-term use, as extra support during sit-to-stand transfers and sitting transfers from the wheelchair to the bed, toilet or car.

DisposableBelt can also be used as support when the user is standing or walking.

DisposableBelt is made of a strong and durable non-woven material and is disposed of after use, when it gets soiled or when the user does not need it anymore, instead of being cleaned. This prevents the spread of multiresistant and other types of contagious bacteria

Art. No.	Product name	Description	Unit
01-6054	DisposableBelt	One size, Medium/ Large	5 Pcs



Turning plates



TurnTable

TurnTable can be used for seated transfers between wheelchair and bed, toilet or chair/armchair. For increased security, the TurnTable can be used together with FlexiBelt/EasyBelt. TurnTable is available in three sizes with or without handle.

Art. No.	Product name	Description	Unit
01-5082	TurnTable	Ø 30 cm/ Ø 11.8", without handle, max 135 kg/297 lbs	1 Pce
01-5083	TurnTable	Ø 38 cm/ Ø 15", without handle, max 135 kg/297 lbs	1 Pce
01-5084	TurnTable	Ø 45 cm/ Ø 17.7", max 250 kg/551 lbs	1 Pce
01-5085	TurnTable	Ø 38 cm/ Ø 15", max 135 kg/297 lbs	1 Pce



TurnSafe

For users with weak legs and poor balance, TurnSafe provides support during transfers between bed, wheelchair, toilet and chair. TurnSafe has a footplate that can be rotated, adjustable leg supports and sturdy handles for both the user and the caregiver. A foot brake on the caregiver's side facilitates braking of the footplate.

Art. No.	Product name	Description	Unit
01-400025544	TurnSafe	Max 150 kg/ 330 lbs	1 Pce

Accessories

01-7215	Padding for leg support	For TurnSafe and ReTurn 7100, 7400, 7500	1 Pair
---------	-------------------------	--	--------



Other support aids



MobiStol

MobiStol is mounted on a chair to make it mobile. The chair can now be used for shorter transfers indoors, such as into and out from a table. For users weighing max 120 kg/ 264 lbs.

Art. No.	Product name	Description	Unit
01-7030	MobiStol	For chairs and armchairs with hard wooden chairlegs of 3 x 3 cm/ 1.2" X 1.2"	1 Set



Support accessories



PantLock

PantLock is fastened around the thigh to prevent the trousers from sliding down too far, for example, during visits to the toilet. This makes it easier to reach down and pull the trousers up afterwards. PantLock has Velcro strips that are passed through a metal eyelet, and then secured to each other.

Art. No.	Product name	Description	Unit
01-8055	PantLock	One size	1 Pce
01-8056	PantLock	One size	10 Pcs





Manual lifting aids



LiftSeat 01-4005

LiftSeat

LiftSeat is developed to be used for manual lifting in a seated position. One common area of application is lifting a user into or out of an aircraft or bus seat. LiftSeat is available in two models, one scoop model and one with divided leg supports for easier application. Both models give the user support under the legs and from the spine base up to the shoulders. And both models are easy to store and simple to take along on trips.

Art. No.	Product name	Description	Unit
01-4005	LiftSeat	Max 120 kg/ 264 lbs	1 Pce
01-4003	LiftSeat, divided leg supports	XXSmall/XSmall, max 150 kg/ 330 lbs	1 Pce
01-4006	LiftSeat, divided leg supports	Small/Medium, max 150 kg/ 330 lbs	1 Pce
01-4007	LiftSeat, divided leg supports	Large/X-Large, max 150 kg/ 330 lbs	1 Pce



LiftSeat 01-4003, 01-4006, 01-4007



FlexiMove 01-4011

FlexiMove

With FlexiMove there is very little need for actual lifting. It is primarily designed for moving higher up in bed, turning and sideways transfers. In these situations, pulling is easier than lifting as the underside is smooth and slides very well in the bed. When lifting can not be avoided, FlexiMove can be used for both sitting and supine transfers. FlexiMove is available in three sizes, all for users weighing max 150 kg/330 lbs.

Art. No.	Product name	Description	Unit
01-4011	FlexiMove	102 x 43 cm/ 40,1" x 16.9", 4 handles	1 Pce
01-4012	FlexiMove	112 x 53 cm/ 44.0" x 20.8", 4 handles	1 Pce
01-4016	FlexiMove	152 x 53 cm/ 59.8" x 20.8", 6 handles	1 Pce



FlexiMove 01-4016



Manual lifting aids



MiniStretcher

MiniStretcher is used for sitting and horizontal transfers, primarily to and from an ambulance stretcher. It is easy to slide in place underneath the patient and has a solid construction with grip-friendly handles to secure safe heavy lifting. The padding and neck support ensures maximum comfort for the patient. The flexible MiniStretcher adapts to the shapes of ambulance stretchers, so they can remain under patients for the entire journey.

Suggested equipment for an ambulance stretcher: MiniStretcher is placed under or on top of the ambulance mattress. The storage pocket for EasyGlide is placed around the frame of the ambulance stretcher under the ambulance mattress. In the storage pocket EasyGlide 01-5050 is stored.

Art. No.	Product name	Description	Unit
01-4017	MiniStretcher	135 x 43 cm/ 53.1" x 16.9" 5 handles, max 200 kg/ 440 lbs	1 Pce



Accessories

01-5050	EasyGlide	75 x 33 cm/ 29.5" x 13.2", creased	1 Pce
01-4014	Storage pocket	For 01-5020, 01-5050	1 Pce



EasyGlide 01-5050



Ceiling lift units



RiseAtlas450M
(sling bar is ordered separately)



Parking Placard
(sling bar is not included)



Cover to charger


RiseAtlas

RiseAtlas is a safe and easy-to-use ceiling lift unit available in several different models and versions with different lifting capacities and with (T) or without (M) built-in transfer motor and many other features and functionalities. When used with the appropriate lifting accessories, RiseAtlas provides secure, comfortable and individually adapted lifting solutions for all situations and user needs. It can be used in all settings and in all types of rooms, both on permanent rail systems and on portable, freestanding lift stands. RiseAtlas offers unusually small measurements, extremely low weight and several unique features as standard. All models of RiseAtlas are supplied with hand control and charger. Sling bar is ordered separately.

Art. No.	Product name	Description	Unit
01-50100042	RiseAtlas450M	SWL 205 kg/450 lbs	1 Pce
01-50100041	RiseAtlas450M	With QuickTrolleySystem, SWL 205 kg/450 lbs	1 Pce
01-50100046	RiseAtlas450M	With switch manoeuvring, SWL 205 kg/450 lbs	1 Pce
01-50100043	RiseAtlas450T	SWL 205 kg/450 lbs	1 Pce
01-50100044	RiseAtlas625M	SWL 285 kg/625 lbs	1 Pce
01-50100047	RiseAtlas625M	With switch manoeuvring, SWL 285 kg/625 lbs	1 Pce
01-50100045	RiseAtlas625T	SWL 285 kg/625 lbs	1 Pce

Accessories

01-70200047	9V battery and cable	For emergency lowering	1 Set
01-70200048	Software for diagnostics	Software	1 Set
01-70200046	QuickTrolley	For RiseAtlas450M, 50100042, 50100041	1 Set
01-50400066	Cover to charger	For RiseAtlas	1 Pce
01-50400048	Parking Placard	For convenient storing of sling bars, manuals and lifting slings	1 Pce
01-70200001	SlingBar S	35 cm / 13.8"	1 Pce
01-70200002	SlingBar M	45 cm / 17.7"	1 Pce
01-70200003	SlingBar L	60 cm / 23.6"	1 Pce
01-70200042	SlingBarSpreader M	Side bars for SlingBar S/M/L	1 Pair



Ceiling lift units



RisePorto450
(trolley is ordered separately)

RisePorto

RisePorto is one of the smallest and lightest portable ceiling lift units available on the market. At the same time, it is both strong and safe, providing secure, comfortable and individually adapted lifting solutions for all situations. RisePorto is available in two models with different lifting capacities.

RisePorto can be used on all types of rails, in all settings and in all types of rooms. The lift unit can easily be moved between different rail systems and rooms, allowing the same lift unit to be used in many locations, giving it a high degree of utilization. RisePorto features a sling bar integrated with the lift unit. RisePorto is very easy to operate, either via the accompanying hand control or directly from the lift unit itself. The trolley(s) is ordered separately.

Art. No.	Product name	Description	Unit
01-50200012	RisePorto300	Complete with hand control and charger, SWL 140 kg / 300 lbs	1 Pce
01-50200011	RisePorto450	Complete with hand control and charger, SWL 205 kg / 450 lbs	1 Pce

Accessories

01-50200008	Trolley RisePorto	For MilkyWay rail system	1 Pce
01-50200004	Room transfer strap	30 cm/ 11.8"	1 Pce
01-70200021	Room transfer strap	70 cm/ 27.6"	1 Pce
01-70200022	Room transfer strap	80 cm/ 31.5"	1 Pce
01-70200023	Room transfer strap	90 cm/ 35.4"	1 Pce
01-70200024	Room transfer strap	100 cm/ 39.4"	1 Pce
01-50200005	Handle for mounting	For RisePorto, 65 cm/ 25.6"	1 Pce
01-70200048	Software for diagnostics	Software	1 Set



Hand control RisePorto

Rail systems



Castor shown with rail and extra locking handles, both of which are sold separately.

Castor and Pollux

Castor is a safe, high-performance and portable freestanding lift stand, which together with a rail and a lift unit, functions as a stationary lift system.

Pollux is a portable free standing lift stand, which together with two rails, BariTurn and two ceiling lift units, functions as a stationary lift system for bariatric users.

Castor and Pollux are made of aluminium, which makes them light in relation to their large loading capacity; Castor 300 kg/660 lbs, Pollux 500 kg/1100 lbs. The design is both compact and discreet while at the same time making Castor and Pollux two of the sturdiest freestanding lift stands on the market. Castor and Pollux are ideal for occasional lifting needs in a specific location, or where it is not possible to install a permanent rail system on the ceiling.

For both Castor and Pollux, the rails are ordered separately, which allows you to choose the rail length and thereby adapt the lift stand's width and the lifting range width to optimally meet lifting needs in different settings. Rails are available in lengths from 2-5 metres.



Pollux shown with rails and BariTurn, both of which are sold separately.

Art. No.	Product name	Description	Unit
01-50300004	Castor	Portable lift stand, without rail, white max 300 kg/ 660 lbs	1 Set
01-50300005	Pollux	Portable lift stand, without rails, white max 500 kg/ 1100 lbs	1 Set
01-50400123	Rail 64-2 W	2 m/ 78"	1 Pce
01-50400157	Rail 64-3 W	3 m/ 118"	1 Pce
01-50400130	Rail 120-4 W	4 m/ 157"	1 Pce
01-50400131	Rail 120-5 W	5 m/ 196"	1 Pce

Accessories

01-50300006	Locking handles	Quick release, for Castor and Pollux	4 Pcs
01-50400180	BariTurn	360° turning trolley for two ceiling lift units, for Pollux	1 Pce



Locking handles



Mobile lifts



MultiLift550

MultiLift550

MultiLift550 is a unique, motor assisted mobile patient lift for all lifting situations and especially useful for heavier users weighing up to 250 kg/550 lbs. MultiLift550 is a truly ergonomic, mobile lifting device requiring no physical exertion from the caregiver. MultiLift550 features powered moving and turning in all directions and at adjustable speed. The large rear drive wheels and the very small turning radius make MultiLift550 very easy to manoeuvre and turn in most spaces and on all types of surfaces including softer floors and carpets. The lift can also easily drive over door thresholds, electric cords and other small obstacles on the floor. MultiLift550 is easily operated via an intuitive and height-adjustable control system consisting of steering handles, a control panel with a display and an infrared remote control. MultiLift550 is supplied complete and ready to use, equipped with a wider sling bar (60 cm / 23.6"), providing a more open lifting sling which is ideal for larger users.

Art. No.	Product name	Description	Unit
01-1011000110	MultiLift550	Complete with remote control, sling bar and charger	1 Pce
Accessories			
01-70200043	SlingBar S, 35 cm / 13.8"	For MultiLift550	1 Pce
01-70200044	SlingBar M, 45 cm / 17.7"	For MultiLift550	1 Pce
01-70200042	SlingBarSpreader M	Side bars for SlingBar S/M/L	1 Pair





Mobile lifts



Eva400EE



Eva450EE



Eva600EE



AmbulationArm

Eva

Eva is a versatile mobile patient lift for all lifting situations. Together with SystemRoMedic's functional and easy-to-use accessories, Eva is the right choice for lifting in connection with transfers to and from the floor, chair or bed, as well as for horizontal lifting and gait training. Eva is available in several models and versions with different lifting capacities and with electrical (EE) or manual (EM) base widening. Eva is also available in special versions with low legs for use with lower beds. All models of Eva are supplied complete with sling bar and hand control. All versions of Eva400 and Eva 450 are equipped with a standard-width sling bar (SlingBar M, 45 cm/ 17.7"). All versions of Eva600 are equipped with a wider sling bar (SlingBar L, 60 cm/ 23.6").

The accessory armrests, AmbulationArm, can be employed by the user during gait training but they can also be used by the carer for manoeuvring the lift.

Art. No.	Product name	Description	Unit
01-60100001	Eva400EE	Max 180 kg/ 400 lbs, aluminum, castor size Ø 100 mm/ 3.9",	1 Pce
01-60100011	Eva400EM	Max 180 kg/ 400 lbs, aluminum, castor size Ø 100 mm/ 3.9",	1 Pce
01-60100002	Eva450EE	Max 205 kg/ 450 lbs, aluminum, castor size Ø 100 mm/ 3.9",	1 Pce
01-60100006	Eva450EE	With low legs Max 205 kg/ 450 lbs, aluminum, castor size front 2 x Ø 35 mm/ 1.4", rear Ø 100 mm/ 3.9"	1 Pce
01-60100007	Eva450EM	With low legs Max 205 kg/ 450 lbs, aluminum, castor size front 2 x Ø 35 mm/ 1.4", rear Ø 100 mm/ 3.9"	1 Pce
01-60100003	Eva600EE	Max 270 kg/ 600 lbs, aluminum, castor size Ø 100 mm/ 3.9"	1 Pce
01-60100010	Eva600EE	With low legs Max 270 kg/ 600 lbs, aluminum, castor size front 2 x Ø 35 mm/ 1.4", rear Ø 100 mm/ 3.9"	1 Pce
Accessories			
01-80100010	AmbulationArm	For Eva mobile lifts	1 Pce
01-70200029	External charger Linak	For Eva450/ Eva600	1 Pce





Rino600EE

Rino600EE is a mobile lift which, together with the right lifting accessories, solves most lifting requirements. It can be used both for seated and horizontal transfers and, thanks to its very high lifting height, it is also excellent for gait training. Rino600EE is made of steel, making it both strong and stable. It can lift users weighing up to 270 kg/600 lbs.

To get more room in a two-point sling bar, SideBar sidebars can be used as an accessory. They open up the sling and provide for a more reclined position.

Art. No.	Product name	Description	Unit
01-60200003	Rino600EE	Max 270 kg/ 600 lbs, steel	1 Pce

Accessories

01-70200013	SideBar	Fits Rino600EE/ Vega505EE	1 Pair
-------------	---------	---------------------------	--------



SideBar




Vega505EE

Vega is a mobile patient lift adapted for safe and easy lifting of users weighing up to 230 kg/ 505 lbs. Vega has a new, unique design with a lift arm that can be both raised to an almost upright position and lowered all the way down to the floor. It is also relatively narrow and compact, with dimensions adapted for use in situations and premises where space is limited. In spite of the handy size, Vega offers an unusually high lifting height and the total lifting range is also extremely large.

Art. No.	Product name	Description	Unit
01-60600002	Vega505EE	Max 230 kg/ 505 lbs, steel, castor size Ø 75 mm/ 3"	1 Pce

Accessories

01-70200013	SideBar	Fits Vega505EE/ Rino600EE	1 Pair
-------------	---------	---------------------------	--------



Lifting accessories



SlingBar

SlingBar is a two-point sling bar in aluminium available in different widths. The sling bars have safety latches which prevent the sling straps from creeping out of the sling bar.

Art. No.	Product name	Description	Unit
01-70200001	SlingBar S	35 cm/ 13.8", max 300 kg/ 660 lbs	1 Pce
01-70200002	SlingBar M	45 cm/ 17.7", max 300 kg/ 660 lbs	1 Pce
01-70200003	SlingBar L	60 cm/ 23.6", max 300 kg/ 660 lbs	1 Pce
01-70200043	SlingBar S	35 cm / 13.8", for MultiLift550	1 Pce
01-70200044	SlingBar M	45 cm / 17.7", for MultiLift550	1 Pce



SlingBarSpreader M,
mounted on SlingBar M.

SlingBarSpreader M

To get more room in a two-point sling bar, SlingBarSpreader M side bars can be used as an accessory. SlingBarSpreader M open up the lifting sling and provide for a more reclined position. SlingBarSpreader M are delivered in pairs and can be used together with all models of SlingBar.

Art. No.	Product name	Description	Unit
01-70200042	SlingBarSpreader M	Side bars for SlingBar S/M/L	1 Pair



SlingBarWrap

SlingBar can be equipped with a soft protection, SlingBarWrap, which is sometimes needed for users with uncontrollable movements. It is attached with Velcro and can easily be removed for cleaning.




SlingBarWrap MiniLift

Art. No.	Product name	Description	Unit
01-70200009	SlingBarWrap	For SlingBar S/M/L	1 Pce
01-70200012	SlingBarWrap MiniLift	For MiniLift200	1 Pce



Lifting accessories




SwiftHook

SwiftHook

With SwiftHook, changing sling bars or connecting a scale becomes very easy.

Art. No.	Product name	Description	Unit
01-70200008	SwiftHook	For SlingBar S/M/L	1 Pce



RoomTransferHook

RoomTransferHook is a quick release to ceiling lift units. Using RoomTransferHook on a two-point sling bar, the user can be transferred between two ceiling lift units mounted on separate rail systems in two separate rooms.


Art. No.	Product name	Description	Unit
01-70200010	RoomTransferHook	Incl. 2 pcs H-adapter + bolts	1 Set

Accessories

01-70200011	H-adapter	Adapter for RoomTransferHook 70200010, 1 Pce for Rise450*/Rise625*
-------------	-----------	--



** Rise450 and Rise625, discontinued*



H-adapter

Lifting accessories



StretcherBar

StretcherBar is a sling bar used together with StretcherSling for horizontal lifting with stationary or mobile lifts. StretcherBar features ten suspension points, five on each side of the user, which provide uniform support under the entire body. The width of the sling is easily adjusted by pulling out or pushing in the suspension hooks.


Art. No.	Product name	Description	Unit
01-70200006	StretcherBar	10 point, max 300 kg/ 660 lbs	1 Pce



Scales

Charder MHS2500 is used together with stationary or mobile lifts for weighing of users.

Art. No.	Product name	Description	Unit
01-70100002	Charder MHS2500	Max 300 kg/ 660 lbs	1 Pce
01-70100003	Charder MHS2500	Max 400 kg/ 880 lbs	1 Pce



BackBoardStraps

BackBoardStraps have been designed for use in combination with the Ferno Scoop Stretcher 065 EXL to enable mechanical lifting with a manual stretcher in emergency care. BackBoardStraps consist of four independent lift straps. There are six alternative lift loops on one side of each strap to facilitate balancing of the stretcher. The BackBoardStraps are looped securely through the handle holes (two on each side), and then connected to the sling bar hooks. SWL 160 kg/352 lbs.

Art. No.	Product name	Description	Unit
01-49482000	BackBoardStraps	For Ferno Scoop Stretcher 065 EXL	4 Pcs



Lifting accessories



BariTurn

BariTurn is a 360° turning trolley for two ceiling lift units. BariTurn is used, in combination with Pollux or a permanently installed MilkyWay rail system, for lifting of extremely large and heavy users, e.g. in connection with sitting and horizontal transfers. BariTurn is also a perfect assistive device for standing and gait training with very heavy users.

Art. No.	Product name	Description	Unit
01-50400180	BariTurn	SWL 500 kg/ 1102 lbs	1 Pce



Consider this when selecting sling size and material

Size

The lifting sling is too small

- if the user's face is too close to the sling bar
- if the leg supports chafe at the groin
- if the lifting sling looks or feels too tight

The lifting sling is too large

- if the user's seat slides out too much in the lower opening
- if the leg supports slide forward towards the hollows of the user's knees

Material

Polyester

- Durable and easy-to-care-for
- Low friction
- Withstands high laundry temperatures

Polyester net

- Flexible and airy
- Water passes through the mesh
- Withstands high laundry temperatures
- Can be left under the user

Plastic net


- Does not absorb water, excellent for lifting in connection with bathing
- Can be disinfected with rubbing alcohol

Disposable material

- Sturdy and durable non-woven
- Instead of laundering, these lifting slings are discarded after use
- Prevents the spread of multi-resistant bacteria

Important details

SystemRoMedic™ offers a wide range of comfortable and functional lifting slings with a safe design and many practical details. Always read the manual before using the lifting sling and don't forget to perform regular function tests.



Lifting
SystemRoMedic™



Lifting slings

Divided leg supports / High back



HighBackSling

Standard model for users who need a lot of support

HighBackSling is designed for use in most common lifting situations, e.g. for sitting transfers between bed and wheelchair, and when the user needs extra head support. HighBackSling can also be used for lifting from a recumbent position on the floor. This standard model suits most users. The high back ensures that the user is supported from the tailbone up over the head and under the legs. HighBackSling is available in several different materials, versions and sizes for users weighing max 300 kg/660 lbs.

Art. No.	Product name	Description	Unit
01-4080002	HighBackSling	Polyester, XX-Small	1 Pce
01-4080003	HighBackSling	Polyester, X-Small	1 Pce
01-4080004	HighBackSling	Polyester, Small	1 Pce
01-4080006	HighBackSling	Polyester, Medium	1 Pce
01-4080007	HighBackSling	Polyester, Large	1 Pce
01-4080008	HighBackSling	Polyester, X-Large	1 Pce
01-4080009	HighBackSling	Polyester, XX-Large	1 Pce
01-4082002	HighBackSling	Polyester net, XX-Small	1 Pce
01-4082003	HighBackSling	Polyester net, X-Small	1 Pce
01-4082004	HighBackSling	Polyester net, Small	1 Pce
01-4082006	HighBackSling	Polyester net, Medium	1 Pce
01-4082007	HighBackSling	Polyester net, Large	1 Pce
01-4082008	HighBackSling	Polyester net, X-Large	1 Pce
01-4082009	HighBackSling	Polyester net, XX-Large	1 Pce
01-4082202	HighBackSling	Polyester net, not reinforced, XX-Small	1 Pce
01-4082203	HighBackSling	Polyester net, not reinforced, X-Small	1 Pce
01-4082204	HighBackSling	Polyester net, not reinforced, Small	1 Pce
01-4082206	HighBackSling	Polyester net, not reinforced, Medium	1 Pce
01-4082207	HighBackSling	Polyester net, not reinforced, Large	1 Pce
01-4082208	HighBackSling	Polyester net, not reinforced, X-Large	1 Pce
01-4082209	HighBackSling	Polyester net, not reinforced, XX-Large	1 Pce
01-4085202	HighBackSling	Plastic net, XX-Small	1 Pce
01-4085203	HighBackSling	Plastic net, X-Small	1 Pce
01-4085204	HighBackSling	Plastic net, Small	1 Pce
01-4085206	HighBackSling	Plastic net, Medium	1 Pce
01-4085207	HighBackSling	Plastic net, Large	1 Pce
01-4085208	HighBackSling	Plastic net, X-Large	1 Pce
01-4085209	HighBackSling	Plastic net, XX-Large	1 Pce



Lifting slings

Divided leg supports / High back



HighBackSling, disposable

Disposable material for optimal hygiene and safety

HighBackSling, disposable, is made of a sturdy, durable non-woven material. Instead of laundering, the lifting sling is discarded after use, when it becomes soiled or when the user no longer needs it. This prevents the spread of multi-resistant and other types of infectious bacteria.

HighBackSling, disposable, is a lifting sling for short-term use and for users weighing up to max 200 kg/440 lbs.

Art. No.	Product name	Description	Unit
0140892004	HighBackSling	Disposable, Small	5 Pcs
01-40892006	HighBackSling	Disposable, Medium	5 Pcs
01-40892007	HighBackSling	Disposable, Large	5 Pcs
01-40892008	HighBackSling	Disposable, X-Large	5 Pcs



ShellHBSling

Unique bowl shape with adjustable head support

ShellHBSling is designed for use in most common lifting situations, e.g. for sitting transfers between bed and wheelchair, and when the user needs extra head support. ShellHBSling can also be used for lifting from a recumbent position on the floor.

The bowl-shaped sling follows the shape of the body and enables an upright seated posture that is easily regulated with the help of alternative lift loops. The adjustable head support gives each user perfect support for the head. The high back ensures that the user is supported from the tailbone up over the head and under the legs. ShellHBSling is available in several sizes, all for users weighing max 300 kg/660 lbs.

Art. No.	Product name	Description	Unit
01-40700002	ShellHBSling	Polyester, XX-Small	1 Pce
01-40700003	ShellHBSling	Polyester, X-Small	1 Pce
01-40700004	ShellHBSling	Polyester, Small	1 Pce
01-40700006	ShellHBSling	Polyester, Medium	1 Pce
01-40700007	ShellHBSling	Polyester, Large	1 Pce
01-40700008	ShellHBSling	Polyester, X-Large	1 Pce
01-40700009	ShellHBSling	Polyester, XX-Large	1 Pce



Lifting slings

Divided leg supports / High back



FlexibleSling

Pliable support for users with deformities

FlexibleSling is a sling that meets the need for a lot of support and, in a flexible way, follows the contours of the body and gives the user a sense of security. FlexibleSling is designed for use in most common lifting situations, e.g. for sitting transfers between bed and wheelchair, and especially for users with deformities and users who need extra head support. FlexibleSling can also be used for lifting from a recumbent position on the floor. The high back ensures that the user is supported from the tailbone up over the head and under the legs. FlexibleSling is available in two versions and several sizes, all for users weighing max 300 kg/660 lbs. FlexibleSling is also available with undivided leg supports and with undivided leg supports and commode opening.

Art. No.	Product name	Description	Unit
01-42520002	FlexibleSling	Polyester net, XX-Small	1 Pce
01-42520003	FlexibleSling	Polyester net, X-Small	1 Pce
01-42520004	FlexibleSling	Polyester net, Small	1 Pce
01-42520006	FlexibleSling	Polyester net, Medium	1 Pce
01-42520007	FlexibleSling	Polyester net, Large	1 Pce
01-42520008	FlexibleSling	Polyester net, X-Large	1 Pce
01-42522002	FlexibleSling	Polyester net, not reinforced, XX-Small	1 Pce
01-42522003	FlexibleSling	Polyester net, not reinforced, X-Small	1 Pce
01-42522004	FlexibleSling	Polyester net, not reinforced, Small	1 Pce
01-42522006	FlexibleSling	Polyester net, not reinforced, Medium	1 Pce
01-42522007	FlexibleSling	Polyester net, not reinforced, Large	1 Pce
01-42522008	FlexibleSling	Polyester net, not reinforced, X-Large	1 Pce



HygieneHBSling

For toileting and for users who need a lot of support

HygieneHBSling is a lifting sling that is easy to apply and makes it easier for the caregiver to dress/undress the user in connection with toileting. HygieneHBSling is designed for users who need a lot of extra support, even for the head. The high back ensures that the user is supported behind the back, up behind the head, around the upper body and under the legs. HygieneHBSling is available in several sizes, all for users weighing max 300 kg/660 lbs.

Art. No.	Product name	Description	Unit
01-43704002	HygieneHBSling	Polyester, XX-Small	1 Pce
01-43704003	HygieneHBSling	Polyester, X-Small	1 Pce
01-43704004	HygieneHBSling	Polyester, Small	1 Pce
01-43704006	HygieneHBSling	Polyester, Medium	1 Pce
01-43704007	HygieneHBSling	Polyester, Large	1 Pce
01-43704008	HygieneHBSling	Polyester, X-Large	1 Pce
01-43704009	HygieneHBSling	Polyester, XX-Large	1 Pce



Lifting slings

Divided leg supports / Shoulder-high back



BasicSling

Standard model for users with good head stability

BasicSling is designed for use in most common lifting situations, e.g. for sitting transfers between bed and wheelchair, and when the patient has poor torso stability but good head stability. This standard model suits most users with good head stability. The shoulder-high back ensures that the user is supported from the tailbone up to the shoulders and under the legs. BasicSling is available in several different materials, versions and sizes, all for users weighing max 300 kg/660 lbs.

Art. No.	Product name	Description	Unit
01-40500002	BasicSling	Polyester, XX-Small	1 Pce
01-40500003	BasicSling	Polyester, X-Small	1 Pce
01-40500004	BasicSling	Polyester, Small	1 Pce
01-40500006	BasicSling	Polyester, Medium	1 Pce
01-40500007	BasicSling	Polyester, Large	1 Pce
01-40500008	BasicSling	Polyester, X-Large	1 Pce
01-40500009	BasicSling	Polyester, XX-Large	1 Pce
01-40520002	BasicSling	Polyester net, XX-Small	1 Pce
01-40520003	BasicSling	Polyester net, X-Small	1 Pce
01-40520004	BasicSling	Polyester net, Small	1 Pce
01-40520006	BasicSling	Polyester net, Medium	1 Pce
01-40520007	BasicSling	Polyester net, Large	1 Pce
01-40520008	BasicSling	Polyester net, X-Large	1 Pce
01-40520009	BasicSling	Polyester net, XX-Large	1 Pce
01-40522002	BasicSling	Polyester net, not reinforced, XX-Small	1 Pce
01-40522003	BasicSling	Polyester net, not reinforced, X-Small	1 Pce
01-40522004	BasicSling	Polyester net, not reinforced, Small	1 Pce
01-40522006	BasicSling	Polyester net, not reinforced, Medium	1 Pce
01-40522007	BasicSling	Polyester net, not reinforced, Large	1 Pce
01-40522008	BasicSling	Polyester net, not reinforced, X-Large	1 Pce
01-40522009	BasicSling	Polyester net, not reinforced, XX-Large	1 Pce
01-40552002	BasicSling	Plastic net, XX-Small	1 Pce
01-40552003	BasicSling	Plastic net, X-Small	1 Pce
01-40552004	BasicSling	Plastic net, Small	1 Pce
01-40552006	BasicSling	Plastic net, Medium	1 Pce
01-40552007	BasicSling	Plastic net, Large	1 Pce
01-40552008	BasicSling	Plastic net, X-Large	1 Pce
01-40552009	BasicSling	Plastic net, XX-Large	1 Pce



Lifting slings

Divided leg supports / Shoulder-high back



AmpSling

For users with high leg amputations

AmpSling is designed for use in most common lifting situations, e.g. for sitting transfers between bed and wheelchair, even when a lifting sling without divided leg supports would have been preferable, for example, when lifting patients with double, high leg amputations. The shoulder-high back ensures that the user is supported from the tailbone up to the shoulders and under the legs.

AmpSling is available in different materials and several sizes, all for users weighing max 300 kg/660 lbs.

Art. No.	Product name	Description	Unit
01-42700003	AmpSling	Polyester, X-Small	1 Pce
01-42700004	AmpSling	Polyester, Small	1 Pce
01-42700006	AmpSling	Polyester, Medium	1 Pce
01-42700007	AmpSling	Polyester, Large	1 Pce
01-42700008	AmpSling	Polyester, X-Large	1 Pce
01-42700009	AmpSling	Polyester, XX-Large	1 Pce
01-42720003	AmpSling	Polyester net, X-Small	1 Pce
01-42720004	AmpSling	Polyester net, Small	1 Pce
01-42720006	AmpSling	Polyester net, Medium	1 Pce
01-42720007	AmpSling	Polyester net, Large	1 Pce
01-42720008	AmpSling	Polyester net, X-Large	1 Pce
01-42720009	AmpSling	Polyester net, XX-Large	1 Pce



A close-up photograph of a blue lifting device, likely a lifting belt or sling. The device is made of a textured, woven fabric. A prominent feature is a wide, bright yellow reflective strip with a black border, positioned horizontally across the lower part of the frame. The device is secured with blue webbing straps and stitching. In the top left corner, a small white tag with the word "RoMedic" is visible. The background is plain white.

Lifting

System **RoMedic**TM

Lifting slings

Divided leg supports / Low back



LowBackSling

Standard model for more active users

LowBackSling is designed for use in most common lifting situations, e.g. for sitting transfers between bed and wheelchair, and when the user has torso stability. This standard model suits most users with torso stability. The low back ensures that the user is supported from the tailbone up to the armpits and under the legs. LowBackSling is available in several materials and sizes, all for users weighing max 300 kg/660 lbs.

Art. No.	Product name	Description	Unit
01-4020002	LowBackSling	Polyester, XX-Small	1 Pce
01-4020003	LowBackSling	Polyester, X-Small	1 Pce
01-4020004	LowBackSling	Polyester, Small	1 Pce
01-4020006	LowBackSling	Polyester, Medium	1 Pce
01-4020007	LowBackSling	Polyester, Large	1 Pce
01-4020008	LowBackSling	Polyester, X-Large	1 Pce
01-4020009	LowBackSling	Polyester, XX-Large	1 Pce
01-4022002	LowBackSling	Polyester net, XX-Small	1 Pce
01-4022003	LowBackSling	Polyester net, X-Small	1 Pce
01-4022004	LowBackSling	Polyester net, Small	1 Pce
01-4022006	LowBackSling	Polyester net, Medium	1 Pce
01-4022007	LowBackSling	Polyester net, Large	1 Pce
01-4022008	LowBackSling	Polyester net, X-Large	1 Pce
01-4022009	LowBackSling	Polyester net, XX-Large	1 Pce
01-40252002	LowBackSling	Plastic net, XX-Small	1 Pce
01-40252003	LowBackSling	Plastic net, X-Small	1 Pce
01-40252004	LowBackSling	Plastic net, Small	1 Pce
01-40252006	LowBackSling	Plastic net, Medium	1 Pce
01-40252007	LowBackSling	Plastic net, Large	1 Pce
01-40252008	LowBackSling	Plastic net, X-Large	1 Pce
01-40252009	LowBackSling	Plastic net, XX-Large	1 Pce



Lifting slings

Divided leg supports / Low back




HygieneSling

For toileting - three versions for more active users

HygieneSling is a lifting sling that is easy to apply and makes it easier for the caregiver to dress/undress the user in connection with toileting. HygieneSling, which is designed for users with adequate muscular tonus. The low back ensures that the user is supported behind the back, up to the armpits and under the legs.

HygieneSling is available in three versions; with or without a belt, or with a belt and extra back support, and in several sizes, all for users weighing max 300 kg/660 lbs.

Art. No.	Product name	Description	Unit
01-43504002	HygieneSling	Polyester, with belt, XX-Small	1 Pce
01-43504003	HygieneSling	Polyester, with belt, X-Small	1 Pce
01-43504004	HygieneSling	Polyester, with belt, Small	1 Pce
01-43504006	HygieneSling	Polyester, with belt, Medium	1 Pce
01-43504007	HygieneSling	Polyester, with belt, Large	1 Pce
01-43504008	HygieneSling	Polyester, with belt, X-Large	1 Pce
01-43504009	HygieneSling	Polyester, with belt, XX-Large	1 Pce
01-43604002	HygieneSling	Polyester, without belt, XX-Small	1 Pce
01-43604003	HygieneSling	Polyester, without belt, X-Small	1 Pce
01-43604004	HygieneSling	Polyester, without belt, Small	1 Pce
01-43604006	HygieneSling	Polyester, without belt, Medium	1 Pce
01-43604007	HygieneSling	Polyester, without belt, Large	1 Pce
01-43604008	HygieneSling	Polyester, without belt, X-Large	1 Pce
01-43604009	HygieneSling	Polyester, without belt, XX-Large	1 Pce
01-43904003	HygieneSling	Polyester, extra back support, X-Small	1 Pce
01-43904004	HygieneSling	Polyester, extra back support, Small	1 Pce
01-43904006	HygieneSling	Polyester, extra back support, Medium	1 Pce
01-43904007	HygieneSling	Polyester, extra back support, Large	1 Pce
01-43904008	HygieneSling	Polyester, extra back support, X-Large	1 Pce



Lifting slings

Divided leg supports / Low back



HygieneLBSling

For toileting - extra support for users with good head and torso stability

HygieneLBSling is a lifting sling that is easy to apply and makes it easier for the caregiver to dress/undress the user in connection with toileting. HygieneLBSling is designed for users with good head and torso stability and adequate muscular tonus. The low back ensures that the user is supported behind the back, up to the armpits, around the upper body and under the legs.

HygieneLBSling is available in several sizes, all for users weighing max 300 kg/660 lbs.

Art. No.	Product name	Description	Unit
01-43804002	HygieneLBSling	Polyester, XX-Small	1 Pce
01-43804003	HygieneLBSling	Polyester, X-Small	1 Pce
01-43804004	HygieneLBSling	Polyester, Small	1 Pce
01-43804006	HygieneLBSling	Polyester, Medium	1 Pce
01-43804007	HygieneLBSling	Polyester, Large	1 Pce
01-43804008	HygieneLBSling	Polyester, X-Large	1 Pce
01-43804009	HygieneLBSling	Polyester, XX-Large	1 Pce



Lifting slings

Undivided leg supports / High back




FlexibleSling, undivided

Gentle and flexible, with maximum support

FlexibleSling, undivided, is a lifting sling that meets the need for a lot of support and, in a flexible way, follows the contours of the body and gives the user a sense of security. FlexibleSling, undivided, is designed for use in most common lifting situations, e.g. for sitting transfers between bed and wheelchair, and especially for users with deformities and users who need extra head support. The high back ensures that the user is supported from the tailbone up over the head and under the legs. FlexibleSling, undivided, is also available with commode opening. This version is designed to facilitate lifting in connection with toileting, especially for users with deformities and who need extra head support. Both versions of FlexibleSling, undivided are available in several sizes, all for users weighing max 300 kg/660 lbs.

FlexibleSling is also available with divided leg supports (page 54).



Art. No.	Product name	Description	Unit
01-42521003	FlexibleSling, undivided	Polyester net, X-Small	1 Pce
01-42521004	FlexibleSling, undivided	Polyester net, Small	1 Pce
01-42521006	FlexibleSling, undivided	Polyester net, Medium	1 Pce
01-42521007	FlexibleSling, undivided	Polyester net, Large	1 Pce
01-42521008	FlexibleSling, undivided	Polyester net, X-Large	1 Pce
01-42521403	FlexibleSling, undivided	Polyester net, with commode opening, X-Small	1 Pce
01-42521404	FlexibleSling, undivided	Polyester net, with commode opening, Small	1 Pce
01-42521406	FlexibleSling, undivided	Polyester net, with commode opening, Medium	1 Pce
01-42521407	FlexibleSling, undivided	Polyester net, with commode opening, Large	1 Pce
01-42521408	FlexibleSling, undivided	Polyester net, with commode opening, X-Large	1 Pce



ClassicHBSling

Gentle for users who need a lot of support

For users who need extra support for the head, ClassicHBSling enables a safe and comfortable seated position and gives the user a sense of security. ClassicHBSling is designed for use in most common lifting situations, e.g. for sitting transfers between bed and wheelchair, and is even suitable for users with high leg amputations. The high back ensures that the user is supported all the way from the hollows of the knees up over the head. ClassicHBSling is available in several sizes, all for users weighing max 300 kg/660 lbs.

Art. No.	Product name	Description	Unit
01-41522003	ClassicHBSling	Polyester net, X-Small	1 Pce
01-41522004	ClassicHBSling	Polyester net, Small	1 Pce
01-41522006	ClassicHBSling	Polyester net, Medium	1 Pce
01-41522007	ClassicHBSling	Polyester net, Large	1 Pce
01-41522008	ClassicHBSling	Polyester net, X-Large	1 Pce



Lifting slings

Undivided leg supports / Shoulder-high back



ClassicSling

Gentle support for users with good head stability

ClassicSling enables a safe and comfortable seated position and gives the user a sense of security. ClassicSling is designed for use in most common lifting situations, e.g. for sitting transfers between bed and wheelchair, and is even suitable for users with high leg amputations. The shoulder-high back ensures that the user is supported from the hollows of the knees up to the shoulders and under the legs. ClassicSling is available in several sizes, all for users weighing max 300 kg/660 lbs.

Art. No.	Product name	Description	Unit
01-41222003	ClassicSling	Polyester net, X-Small	1 Pce
01-41222004	ClassicSling	Polyester net, Small	1 Pce
01-41222006	ClassicSling	Polyester net, Medium	1 Pce
01-41222007	ClassicSling	Polyester net, Large	1 Pce
01-41222008	ClassicSling	Polyester net, X-Large	1 Pce



BariSling

For extremely large and heavy users - over 200 kg/440 lbs


BariSling is a lifting sling that is designed for lifting users weighing 200 kg/440lbs or more. It can be used in combination with both mobile and stationary lifts, with a four-point sling bar or with BariTurn. BariSling has six suspension points; two loops at the back and four below, under the legs, to distribute pressure safely and comfortably. The shoulder-high back ensures that the user is supported from the hollows of the knees up to the shoulders and under the legs. BariSling is available in one size for users weighing max 570 kg/1250 lbs.

Art. No.	Product name	Description	Unit
01-47001010	BariSling	Polyester, XXX-Large	1 Pce



Lifting slings


Horizontal lifting / Turning



StretcherSling

StretcherSling is a lifting sling that is used in combination with mobile lifts or stationary lift systems and SystemRoMedic's StretcherBar, a specially designed sling bar for horizontal lifting and transferring users in a lying position. StretcherSling has ten suspension points; five on either side of the user, giving the user evenly distributed support under the entire body. StretcherSling is available in different materials and one size, all for users weighing up to 300 kg/660 lbs.

Art. No.	Product name	Description	Unit
01-46502007	StretcherSling	Polyester, one size	1 Pce
01-46592007	StretcherSling	Disposable, one size	5 Pcs



SafeHandlingSheet

SafeHandlingSheet is a sheet that is used for horizontal lifting or for positioning a user in bed without subjecting the caregiver to unnecessary strain. SafeHandlingSheet is used in combination with mobile lifts or stationary lift systems and a two-point sling bar. SafeHandlingSheet has sixteen lift straps, eight on either side of the user, with three alternative lift loops each, which enables optimal adaptation to the user and transfer situation. SafeHandlingSheet is available in several materials and sizes.

Art. No.	Product name	Description	Unit
01-47522007	SafeHandlingSheet	Polyester net, Large, max 300 kg/660 lbs	1 Pce
01-47542007	SafeHandlingSheet	Polyester net/WendyLett4Way, Large, max 300 kg/660 lbs	1 Pce
01-47502008	SafeHandlingSheet	Polyester, X-Large, max 400 kg/900 lbs	1 Pce
01-47502009	SafeHandlingSheet	Polyester, reinforced, XX-Large, max 570 kg/1250 lbs	1 Pce
01-47592006	SafeHandlingSheet	Disposable, Medium, max 300 kg/660 lbs	3 Pcs



TurningSling

TurningSling is a lifting sling for short-term use. TurningSling is used, in combination with a two-point sling bar mounted on a stationary or mobile lift, for turning users in bed. TurningSling also facilitates hygiene, treatment of sores and dressing. TurningSling is available in one size, for users weighing up to 250 kg/ 550 lbs.

Art. No.	Product name	Description	Unit
01-49592000	TurningSling	Disposable, one size	5 Pcs



Lifting slings

Rising / Standing / Walking



ThoraxSling

ThoraxSling

ThoraxSling is a specially adapted sit-to-stand vest that is used in combination with SystemRo-Medic's sit-to-stand lifts MiniLift200 and MiniLift125. It is easy for the caregiver to apply, while giving the user an added sense of security during rising and gait training. ThoraxSling also facilitates toileting. ThoraxSling is available in two models and in several sizes, all for users weighing max 300 kg/660 lbs.

ThoraxSling, with seat support is also available in a disposable material. This version of ThoraxSling is made of a sturdy, durable non-woven material. Instead of laundering, the lifting vest is discarded after use, when it becomes soiled or when the user no longer needs it. This prevents the spread of multi-resistant and other types of infectious bacteria. ThoraxSling, disposable, is a lifting vest for short-term use and for users weighing up to max 200 kg/440 lbs.



ThoraxSling,
with seat support



ThoraxSling,
with seat support,
disposable

Art. No.	Product name	Description	Unit
01-45500003	ThoraxSling	Polyester, X-Small	1 Pce
01-45500004	ThoraxSling	Polyester, Small	1 Pce
01-45500006	ThoraxSling	Polyester, Medium	1 Pce
01-45500007	ThoraxSling	Polyester, Large	1 Pce
01-45500008	ThoraxSling	Polyester, X-Large	1 Pce
01-45500009	ThoraxSling	Polyester, XX-Large	1 Pce
01-45600004	ThoraxSling	With seat support, polyester, Small	1 Pce
01-45600006	ThoraxSling	With seat support, polyester, Medium	1 Pce
01-45600007	ThoraxSling	With seat support, polyester, Large	1 Pce
01-45600008	ThoraxSling	With seat support, polyester, X-Large	1 Pce
01-45690004	ThoraxSling	With seat support, disposable, Small	5 Pcs
01-45690006	ThoraxSling	With seat support, disposable, Medium	5 Pcs
01-45690007	ThoraxSling	With seat support, disposable, Large	5 Pcs
01-45690008	ThoraxSling	With seat support, disposable, X-Large	5 Pcs



Lifting slings

Rising / Standing / Walking



StandUpSling

StandUpSling is a sit-to-stand vest that is easy for the caregiver to apply, while giving the user an added sense of security. It is designed to provide support behind the back during rising with a sit-to-stand lift such as SystemRoMedic's RiseUp600EE, and to facilitate toileting. StandUpSling enables a comfortable sit-to-stand experience, despite reduced function. The user should have a degree of torso stability.

StandUpSling is available in several sizes, all for users weighing up to 300 kg/660 lbs.

For pressure-sensitive users, StandUpSling Padding is available as an accessory. StandUpSling Padding consists of separate padded cases that give the lift vest an extra layer of soft material and are easy for the caregiver to attach to the StandUpSling.

Art. No.	Product name	Description	Unit
01-45200004	StandUpSling	Polyester, Small	1 Pce
01-45200006	StandUpSling	Polyester, Medium	1 Pce
01-45200007	StandUpSling	Polyester, Large	1 Pce
01-45200008	StandUpSling	Polyester, X-Large	1 Pce
01-45200009	StandUpSling	Polyester, XX-Large	1 Pce

Accessories

01-45273306	StandUpSling Padding	Small/ Medium	1 Pce
01-45273307	StandUpSling Padding	Large/ X-Large	1 Pce



Lifting slings

Rising / Standing / Walking



StandingVest

StandingVest is a lift vest that is specially designed for lifting in connection with standing training when the user uses a standing shell. StandingVest provides support in the groin, around the upper body and all the way up behind the head. It is held in place by a leg harness at the groin.

StandingVest is available in different sizes, all for users weighing max 300 kg/660 lbs.

Art. No.	Product name	Description	Unit
01-44503003	StandingVest	Polyester, X-Small	1 Pce
01-44503004	StandingVest	Polyester, Small	1 Pce
01-44503006	StandingVest	Polyester, Medium	1 Pce
01-44503007	StandingVest	Polyester, Large	1 Pce
01-44503008	StandingVest	Polyester, X-Large	1 Pce



WalkingVest

WalkingVest is a vest designed for rising and gait training with the help of a stationary or mobile lift. WalkingVest provides support at the groin and is kept in place with a padded leg harness.

WalkingVest is available in different sizes, all for users weighing max 300 kg/660 lbs.

Art. No.	Product name	Description	Unit
01-44200002	WalkingVest	Polyester, XX-Small	1 Pce
01-44200003	WalkingVest	Polyester, X-Small	1 Pce
01-44200004	WalkingVest	Polyester, Small	1 Pce
01-44200006	WalkingVest	Polyester, Medium	1 Pce
01-44200007	WalkingVest	Polyester, Large	1 Pce
01-44200008	WalkingVest	Polyester, X-Large	1 Pce
01-44200009	WalkingVest	Polyester, XX-Large	1 Pce



Lifting slings

Rising / Standing / Walking



BariVest

For activating extremely large and heavy users

BariVest is a vest that is designed for rising, standing training and gait training with users weighing 180 kg/397 lbs or more. It can be used in combination with both mobile and stationary lifts, with a four-point sling bar or with BariTurn. BariVest consists of a back section, a front section and a narrower section in the middle to enable a comfortable fit in the groin. It has two sling bar loops, both front and back, as well as belts at three levels.

BariVest is available in different sizes, all for users weighing max 570 kg/1250 lbs.

Art. No.	Product name	Description	Unit
01-47103010	BariVest	Polyester, XXX-Large	1 Pce
01-47103011	BariVest	Polyester, XXXX-Large	1 Pce

Accessories

01-49692000	BariVest cover	Cover for the middle part of BariVest, disposable, one size	5 Pcs
-------------	----------------	---	-------



Lifting slings

Lift straps



LimbSling

LimbSling

LimbSling is designed to comfortably and conveniently assist lifting of a user's leg or arm, for example, to gain better access for treatment and dressing sores. LimbSling is very practical for large and heavy users. LimbSling is applied by placing it under the limb to be lifted and connecting the strap loops to the sling bar.


LimbSling is available in one size and two different materials.

Art. No.	Product name	Description	Unit
01-49200006	LimbSling	Polyester, max 300 kg/ 660 lbs	1 Pce
01-49292000	LimbSling	Disposable, max 150 kg/ 330 lbs	5 Pcs



LimbSling,
disposable

Size guide for lifting slings



Divided leg supports/High back



Art. No.	Product name	Material / Description	Size	Rec user weight (kg)	Sling seat measurement (cm)	Sling overall height (cm)
01-40800002	HighBackSling	polyester	XXS	10-25 / 22-55 lbs	20 / 8"	77 / 30"
01-40800003	HighBackSling	polyester	XS	20-35 / 44-77 lbs	25 / 10"	85 / 33"
01-40800004	HighBackSling	polyester	S	30-55 / 66-122 lbs	35 / 14"	91 / 36"
01-40800006	HighBackSling	polyester	M	50-90 / 110-198 lbs	45 / 18"	94 / 37"
01-40800007	HighBackSling	polyester	L	80-110 / 176-242 lbs	56 / 22"	98 / 39"
01-40800008	HighBackSling	polyester	XL	110-160 / 242-353 lbs	62 / 24"	98 / 39"
01-40800009	HighBackSling	polyester	XXL	130-200 / 287-441 lbs	70 / 28"	102 / 40"
01-40820002	HighBackSling	polyester net	XXS	10-25 / 22-55 lbs	20 / 8"	77 / 30"
01-40820003	HighBackSling	polyester net	XS	20-35 / 44-77 lbs	25 / 10"	85 / 33"
01-40820004	HighBackSling	polyester net	S	30-55 / 66-122 lbs	35 / 14"	91 / 36"
01-40820006	HighBackSling	polyester net	M	50-90 / 110-198 lbs	45 / 18"	94 / 37"
01-40820007	HighBackSling	polyester net	L	80-110 / 176-242 lbs	56 / 22"	98 / 39"
01-40820008	HighBackSling	polyester net	XL	110-160 / 242-353 lbs	62 / 24"	98 / 39"
01-40820009	HighBackSling	polyester net	XXL	130-200 / 287-441 lbs	70 / 28"	102 / 40"
01-40822002	HighBackSling	polyester net, not reinf.	XXS	10-25 / 22-55 lbs	20 / 8"	77 / 30"
01-40822003	HighBackSling	polyester net, not reinf.	XS	20-35 / 44-77 lbs	25 / 10"	85 / 33"
01-40822004	HighBackSling	polyester net, not reinf.	S	30-55 / 66-122 lbs	35 / 14"	91 / 36"
01-40822006	HighBackSling	polyester net, not reinf.	M	50-90 / 110-198 lbs	45 / 18"	94 / 37"
01-40822007	HighBackSling	polyester net, not reinf.	L	80-110 / 176-243 lbs	56 / 22"	98 / 39"
01-40822008	HighBackSling	polyester net, not reinf.	XL	110-160 / 242-353 lbs	62 / 24"	98 / 39"
01-40822009	HighBackSling	polyester net, not reinf.	XXL	130-200 / 287-441 lbs	70 / 28"	102 / 40"
01-40852002	HighBackSling	plastic net	XXS	10-25 / 22-55 lbs	20 / 8"	77 / 30"
01-40852003	HighBackSling	plastic net	XS	20-35 / 44-77 lbs	25 / 10"	85 / 33"
01-40852004	HighBackSling	plastic net	S	30-55 / 66-122 lbs	35 / 14"	91 / 36"
01-40852006	HighBackSling	plastic net	M	50-90 / 110-198 lbs	45 / 18"	94 / 37"
01-40852007	HighBackSling	plastic net	L	80-110 / 176-242 lbs	56 / 22"	98 / 39"
01-40852008	HighBackSling	plastic net	XL	110-160 / 242-353 lbs	62 / 24"	98 / 39"
01-40852009	HighBackSling	plastic net	XXL	130-200 / 287-441 lbs	70 / 28"	102 / 40"
01-40892004	HighBackSling	disposable	S	30-55 / 66-122 lbs	35 / 14"	91 / 36"
01-40892006	HighBackSling	disposable	M	50-90 / 110-198 lbs	45 / 18"	94 / 37"
01-40892007	HighBackSling	disposable	L	80-110 / 176-242 lbs	56 / 22"	98 / 39"
01-40892008	HighBackSling	disposable	XL	110-160 / 242-353 lbs	62 / 24"	98 / 39"
01-40700002	ShellHBSling	polyester	XXS	10-30 / 22-66 lbs	17 / 7"	80 / 31"
01-40700003	ShellHBSling	polyester	XS	20-35 / 44-77 lbs	25 / 10"	85 / 33"
01-40700004	ShellHBSling	polyester	S	30-55 / 66-122 lbs	35 / 14"	94 / 37"
01-40700006	ShellHBSling	polyester	M	50-90 / 110-198 lbs	40 / 16"	97 / 38"
01-40700007	ShellHBSling	polyester	L	80-120 / 176-265 lbs	45 / 18"	98 / 39"
01-40700008	ShellHBSling	polyester	XL	110-160 / 242-353 lbs	55 / 22"	100 / 39"
01-40700009	ShellHBSling	polyester	XXL	150-210 / 331-463 lbs	65 / 26"	100 / 39"
01-42520002	FlexibleSling	polyester net	XXS	10-30 / 22-66 lbs	20 / 8"	75 / 30"
01-42520003	FlexibleSling	polyester net	XS	20-35 / 44-77 lbs	25 / 10"	81 / 32"
01-42520004	FlexibleSling	polyester net, not reinf.	S	25-55 / 55-121 lbs	35 / 14"	90 / 35"
01-42520006	FlexibleSling	polyester net	M	50-75 / 110-165 lbs	45 / 18"	102 / 40"
01-42520007	FlexibleSling	polyester net	L	70-110 / 154-242 lbs	55 / 22"	102 / 40"
01-42520008	FlexibleSling	polyester net	XL	100-130 / 220-287 lbs	65 / 26"	104 / 41"
01-42522002	FlexibleSling	polyester net, not reinf.	XXS	10-30 / 22-66 lbs	20 / 8"	75 / 30"
01-42522003	FlexibleSling	polyester net, not reinf.	XS	20-35 / 44-77 lbs	25 / 10"	81 / 32"
01-42522004	FlexibleSling	polyester net, not reinf.	S	25-55 / 55-121 lbs	35 / 14"	90 / 35"
01-42522006	FlexibleSling	polyester net, not reinf.	M	50-75 / 110-165 lbs	45 / 18"	102 / 40"
01-42522007	FlexibleSling	polyester net, not reinf.	L	70-110 / 154-242 lbs	55 / 22"	102 / 40"
01-42522008	FlexibleSling	polyester net, not reinf.	XL	100-130 / 220-287 lbs	65 / 26"	104 / 41"



PLEASE NOTE!

To determine the most appropriate model, size and material for the individual user, trial fitting of lifting slings is always essential. The size guide's measurements and recommendations are only intended as an aid to selecting lifting slings. The lifting sling's seat width gives an indication of how wide the user's hips (trochanter width) can be. The user's trochanter width must be within +/- 10 cm/3.9" of the lifting sling's seat width. The lift vest's width gives an indication of the girth of the user's trunk. The user's trunk girth must be within +/- 10 cm/3.9" of the lift vest's width. The lifting sling is too large if the user's seat slides out too far at the lower opening, or if the leg supports slide forward towards the hollows of the user's knees. The lifting sling is too small if the user's face is too close to the slingbar, if the leg supports chafe at the groin, or if the sling appears to be too tight or feels too tight.

Divided leg supports/Shoulder-high back



Art. No.	Product name	Material / Description	Size	Rec user weight (kg)	Sling seat measurement (cm)	Sling overall height (cm)
01-40500002	BasicSling	polyester	XXS	10-25 / 22-55 lbs	21 / 8"	63 / 25"
01-40500003	BasicSling	polyester	XS	20-35 / 44-77 lbs	27 / 11"	65 / 26"
01-40500004	BasicSling	polyester	S	35-55 / 77-121 lbs	37 / 15"	77 / 30"
01-40500006	BasicSling	polyester	M	45-80 / 99-176 lbs	50 / 20"	89 / 35"
01-40500007	BasicSling	polyester	L	75-110 / 165-242 lbs	60 / 24"	95 / 37"
01-40500008	BasicSling	polyester	XL	100-150 / 220-330 lbs	65 / 26"	98 / 39"
01-40500009	BasicSling	polyester	XXL	140-200 / 309-441 lbs	75 / 30"	98 / 39"
01-40520002	BasicSling	polyester net	XXS	10-25 / 22-55 lbs	21 / 8"	63 / 25"
01-40520003	BasicSling	polyester net	XS	20-35 / 44-77 lbs	27 / 11"	65 / 26"
01-40520004	BasicSling	polyester net	S	35-55 / 77-121 lbs	37 / 15"	77 / 30"
01-40520006	BasicSling	polyester net	M	45-80 / 99-176 lbs	50 / 20"	89 / 35"
01-40520007	BasicSling	polyester net	L	75-110 / 165-242 lbs	60 / 24"	95 / 37"
01-40520008	BasicSling	polyester net	XL	100-150 / 220-330 lbs	65 / 26"	98 / 39"
01-40520009	BasicSling	polyester net	XXL	140-200 / 309-441 lbs	75 / 30"	98 / 39"
01-40522002	BasicSling	polyester net, not reinf.	XXS	10-25 / 22-55 lbs	21 / 8"	63 / 25"
01-40522003	BasicSling	polyester net, not reinf.	XS	20-35 / 44-77 lbs	27 / 11"	65 / 26"
01-40522004	BasicSling	polyester net, not reinf.	S	35-55 / 77-121 lbs	37 / 15"	77 / 30"
01-40522006	BasicSling	polyester net, not reinf.	M	45-80 / 99-176 lbs	50 / 20"	89 / 35"
01-40522007	BasicSling	polyester net, not reinf.	L	75-110 / 165-242 lbs	60 / 24"	95 / 37"
01-40522008	BasicSling	polyester net, not reinf.	XL	100-150 / 220-330 lbs	65 / 26"	98 / 39"
01-40522009	BasicSling	polyester net, not reinf.	XXL	140-200 / 309-441 lbs	75 / 30"	98 / 39"
01-40552002	BasicSling	plastic net	XXS	10-25 / 22-55 lbs	21 / 8"	63 / 25"
01-40552003	BasicSling	plastic net	XS	20-35 / 44-77 lbs	27 / 11"	65 / 26"
01-40552004	BasicSling	plastic net	S	35-55 / 77-121 lbs	37 / 15"	77 / 30"
01-40552006	BasicSling	plastic net	M	45-80 / 99-176 lbs	50 / 20"	89 / 35"
01-40552007	BasicSling	plastic net	L	75-110 / 165-242 lbs	60 / 24"	95 / 37"
01-40552008	BasicSling	plastic net	XL	100-150 / 220-330 lbs	65 / 26"	98 / 39"
01-40552009	BasicSling	plastic net	XXL	140-200 / 309-441 lbs	75 / 30"	98 / 39"
01-42700003	AmpSling	polyester	XS	20-40 / 44-88 lbs	27 / 11"	85 / 33"
01-42700004	AmpSling	polyester	S	35-60 / 77-133 lbs	37 / 15"	90 / 35"
01-42700006	AmpSling	polyester	M	50-80 / 110-176 lbs	45 / 18"	97 / 38"
01-42700007	AmpSling	polyester	L	70-110 / 154-242 lbs	50 / 20"	100 / 39"
01-42700008	AmpSling	polyester	XL	100-130 / 220-287 lbs	60 / 24"	100 / 39"
01-42700009	AmpSling	polyester	XXL	120-180 / 265-397 lbs	70 / 28"	100 / 39"
01-42720003	AmpSling	polyester net	XS	20-40 / 44-88 lbs	27 / 11"	85 / 33"
01-42720004	AmpSling	polyester net	S	35-60 / 77-133 lbs	37 / 15"	90 / 35"
01-42720006	AmpSling	polyester net	M	50-80 / 110-176 lbs	45 / 18"	97 / 38"
01-42720007	AmpSling	polyester net	L	70-110 / 154-242 lbs	50 / 20"	100 / 39"
01-42720008	AmpSling	polyester net	XL	100-130 / 220-287 lbs	60 / 24"	100 / 39"
01-42720009	AmpSling	polyester net	XXL	120-180 / 265-397 lbs	70 / 28"	100 / 39"



Divided leg supports/Low back



Art. No.	Product name	Material	Size	Rec user weight (kg)	Sling seat measurement (cm)	Sling overall height (cm)
01-40200002	LowBackSling	polyester	XXS	10-25 / 22-55 lbs	21 / 8"	60 / 24"
01-40200003	LowBackSling	polyester	XS	20-35 / 44-77 lbs	27 / 11"	60 / 24"
01-40200004	LowBackSling	polyester	S	30-50 / 66-110 lbs	37 / 15"	77 / 30"
01-40200006	LowBackSling	polyester	M	45-80 / 99-176 lbs	47 / 18"	85 / 33"
01-40200007	LowBackSling	polyester	L	75-110 / 165-242 lbs	57 / 22"	90 / 35"
01-40200008	LowBackSling	polyester	XL	100-140 / 220-309 lbs	67 / 26"	92 / 36"
01-40200009	LowBackSling	polyester	XXL	130-190 / 287-419 lbs	77 / 30"	92 / 36"
01-40220002	LowBackSling	polyester net	XXS	10-25 / 22-55 lbs	21 / 8"	60 / 24"
01-40220003	LowBackSling	polyester net	XS	20-35 / 44-77 lbs	27 / 11"	60 / 24"
01-40220004	LowBackSling	polyester net	S	30-50 / 66-110 lbs	37 / 15"	77 / 30"
01-40220006	LowBackSling	polyester net	M	45-80 / 99-176 lbs	47 / 18"	85 / 33"
01-40220007	LowBackSling	polyester net	L	75-110 / 165-242 lbs	57 / 22"	90 / 35"
01-40220008	LowBackSling	polyester net	XL	100-140 / 220-309 lbs	67 / 26"	92 / 36"
01-40220009	LowBackSling	polyester net	XXL	130-190 / 287-419 lbs	77 / 30"	92 / 36"
01-40252002	LowBackSling	plastic net	XXS	10-25 / 22-55 lbs	21 / 8"	60 / 24"
01-40252003	LowBackSling	plastic net	XS	20-35 / 44-77 lbs	27 / 11"	60 / 24"
01-40252004	LowBackSling	plastic net	S	30-50 / 66-110 lbs	37 / 15"	77 / 30"
01-40252006	LowBackSling	plastic net	M	45-80 / 99-176 lbs	47 / 18"	85 / 33"
01-40252007	LowBackSling	plastic net	L	75-110 / 165-242 lbs	57 / 22"	90 / 35"
01-40252008	LowBackSling	plastic net	XL	100-140 / 220-309 lbs	67 / 26"	92 / 36"
01-40252009	LowBackSling	plastic net	XXL	130-190 / 287-419 lbs	77 / 30"	92 / 36"

Undivided leg supports/High and shoulder-high back



Art. No.	Product name	Material	Size	Rec user weight (kg)	Sling seat measurement (cm)	Sling overall height (cm)
01-42521003	FlexibleSling, undivided	polyester net	XS	20-40 / 44-88 lbs	75 / 30"	80 / 31"
01-42521004	FlexibleSling, undivided	polyester net	S	30-70 / 66-154 lbs	85 / 33"	85 / 33"
01-42521006	FlexibleSling, undivided	polyester net	M	60-90 / 132-198 lbs	94 / 37"	93 / 37"
01-42521007	FlexibleSling, undivided	polyester net	L	80-130 / 176-287 lbs	105 / 41"	96 / 38"
01-42521008	FlexibleSling, undivided	polyester net	XL	110-180 / 242-397 lbs	115 / 45"	98 / 39"
01-41522003	ClassicHBSling	polyester net	XS	20-40 / 44-88 lbs	30 / 12"	95 / 37"
01-41522004	ClassicHBSling	polyester net	S	30-60 / 66-132 lbs	40 / 16"	98 / 39"
01-41522006	ClassicHBSling	polyester net	M	55-85 / 121-187 lbs	48 / 19"	106 / 42"
01-41522007	ClassicHBSling	polyester net	L	80-120 / 176-265 lbs	65 / 26"	110 / 43"
01-41522008	ClassicHBSling	polyester net	XL	110-140 / 242-309 lbs	80 / 31"	115 / 45"
01-41222003	ClassicSling	polyester net	XS	20-40 / 44-88 lbs	30 / 12"	80 / 31"
01-41222004	ClassicSling	polyester net	S	30-60 / 66-132 lbs	38 / 15"	90 / 35"
01-41222006	ClassicSling	polyester net	M	55-85 / 121-187 lbs	46 / 18"	95 / 37"
01-41222007	ClassicSling	polyester net	L	80-120 / 176-265 lbs	63 / 25"	98 / 38"
01-41222008	ClassicSling	polyester net	XL	110-140 / 242-309 lbs	80 / 31"	100 / 39"
01-47001010	BariSling	polyester net	XXXL	200-570 / 441-1257 lbs	-	-

To and from toilet/High back



Art. No.	Product name	Material	Size	Rec user weight (kg)	Sling seat measurement (cm)	Sling overall height (cm)
01-43704002	HygieneHBSling	polyester	XXS	10-30 / 22-66 lbs	25 / 10"	92 / 36"
01-43704003	HygieneHBSling	polyester	XS	25-45 / 55-99 lbs	28 / 11"	92 / 36"
01-43704004	HygieneHBSling	polyester	S	40-60 / 88-132 lbs	30 / 12"	96 / 38"
01-43704006	HygieneHBSling	polyester	M	55-80 / 121-176 lbs	35 / 14"	102 / 40"
01-43704007	HygieneHBSling	polyester	L	75-110 / 165-242 lbs	41 / 16"	102 / 40"
01-43704008	HygieneHBSling	polyester	XL	90-130 / 198-287 lbs	46 / 18"	104 / 41"
43704009	HygieneHBSling	polyester	XXL	120-170 / 265-376 lbs	60 / 24"	104 / 41"
01-42521403	FlexibleSling, undivided*	polyester net	XS	20-40 / 44-88 lbs	75 / 30"	80 / 31"
01-42521404	FlexibleSling, undivided*	polyester net	S	30-70 / 66-154 lbs	85 / 33"	85 / 33"
01-42521406	FlexibleSling, undivided*	polyester net	M	60-90 / 132-198 lbs	94 / 37"	93 / 37"
01-42521407	FlexibleSling, undivided*	polyester net	L	80-130 / 176-287 lbs	105 / 41"	96 / 38"
01-42521408	FlexibleSling, undivided*	polyester net	XL	110-180 / 242-397 lbs	115 / 45"	98 / 38"

*) with commode opening

To and from toilet/Low back



Art. No.	Product name	Material / Description	Size	Rec user weight (kg)	Sling seat measurement (cm)	Sling overall height (cm)
01-43504002	HygieneSling	polyester, with belt	XXS	10-30 / 22-66 lbs	25 / 10"	50 / 20"
01-43504003	HygieneSling	polyester, with belt	XS	25-45 / 55-99 lbs	35 / 14"	65 / 26"
01-43504004	HygieneSling	polyester, with belt	S	40-60 / 88-132 lbs	48 / 19"	75 / 30"
01-43504006	HygieneSling	polyester, with belt	M	55-85 / 121-187 lbs	55 / 22"	87 / 35"
01-43504007	HygieneSling	polyester, with belt	L	75-100 / 165-220 lbs	60 / 24"	93 / 37"
01-43504008	HygieneSling	polyester, with belt	XL	90-130 / 198-287 lbs	80 / 31"	95 / 37"
01-43504009	HygieneSling	polyester, with belt	XXL	120-170 / 265-376 lbs	90 / 35"	95 / 37"
01-43604002	HygieneSling	polyester, without belt	XXS	10-30 / 22-66 lbs	25 / 10"	50 / 20"
01-43604003	HygieneSling	polyester, without belt	XS	25-45 / 55-99 lbs	35 / 14"	65 / 26"
01-43604004	HygieneSling	polyester, without belt	S	40-60 / 88-132 lbs	48 / 19"	75 / 30"
01-43604006	HygieneSling	polyester, without belt	M	55-85 / 121-187 lbs	55 / 22"	87 / 35"
01-43604007	HygieneSling	polyester, without belt	L	75-100 / 165-220 lbs	60 / 24"	93 / 37"
01-43604008	HygieneSling	polyester, without belt	XL	90-130 / 198-287 lbs	80 / 31"	95 / 37"
01-43604009	HygieneSling	polyester, without belt	XXL	120-170 / 265-376 lbs	90 / 35"	95 / 37"
01-43904003	HygieneSling	polyester, extra back support	XS	25-45 / 55-99 lbs	35 / 14"	65 / 26"
01-43904004	HygieneSling	polyester, extra back support	S	40-60 / 88-132 lbs	48 / 19"	75 / 30"
01-43904006	HygieneSling	polyester, extra back support	M	55-85 / 121-187 lbs	55 / 22"	87 / 35"
01-43904007	HygieneSling	polyester, extra back support	L	75-100 / 165-220 lbs	60 / 24"	93 / 37"
01-43904008	HygieneSling	polyester, extra back support	XL	90-130 / 198-287 lbs	80 / 31"	95 / 37"
01-43804002	HygieneLBSling	polyester	XXS	10-30 / 22-66 lbs	25 / 10"	92 / 36"
01-43804003	HygieneLBSling	polyester	XS	25-45 / 55-99 lbs	28 / 11"	92 / 36"
01-43804004	HygieneLBSling	polyester	S	40-60 / 88-132 lbs	30 / 12"	96 / 38"
01-43804006	HygieneLBSling	polyester	M	55-80 / 121-176 lbs	35 / 14"	102 / 40"
01-43804007	HygieneLBSling	polyester	L	75-110 / 165-242 lbs	41 / 16"	102 / 40"
01-43804008	HygieneLBSling	polyester	XL	90-130 / 198-287 lbs	46 / 18"	104 / 41"
01-43804009	HygieneLBSling	polyester	XXL	120-170 / 265-376 lbs	60 / 24"	104 / 41"

Horizontal lifting/Turning



Art. No.	Product name	Material / Description	Size	Rec user weight (kg)	Lenght (cm)	Width (cm)
01-46502007	StretcherSling	polyester	one size	50-200 / 110-441 lbs	180 / 71"	86 / 34"
01-46592007	StretcherSling	disposable	one size	50-200 / 110-441 lbs	180 / 71"	86 / 34"
01-47522007	SafeHandlingSheet	polyester net	L	50-300 / 110-661 lbs	240 / 94"	115 / 45"
01-47542007	SafeHandlingSheet	pol.net/WL4W	L	50-300 / 110-661 lbs	240 / 94"	115 / 45"
01-47502008	SafeHandlingSheet	polyester	XL	50-400 / 110-881 lbs	240 / 94"	115 / 45"
01-47502009	SafeHandlingSheet	polyester, reinf.	XXL	50-570 / 110-1257 lbs	240 / 94"	130 / 51"
01-47592006	SafeHandlingSheet	disposable	M	50-300 / 110-661 lbs	200 / 79"	120 / 47"
01-49592000	TurningSling	disposable	one size	50-250 / 110-551 lbs		

Rising/Standing with MiniLift



Art. No.	Product name	Material / Description	Size	Rec user weight (kg)	Lift vest width (cm)
01-45500003	ThoraxSling	polyester	XS	20-50 / 44-110 lbs	85 / 33"
01-45500004	ThoraxSling	polyester	S	30-60 / 66-132 lbs	94 / 37"
01-45500006	ThoraxSling	polyester	M	50-80 / 110-176 lbs	104 / 41"
01-45500007	ThoraxSling	polyester	L	70-110 / 154-242 lbs	124 / 49"
01-45500008	ThoraxSling	polyester	XL	100-150 / 220-330 lbs	146 / 57"
01-45500009	ThoraxSling	polyester	XXL	140-200 / 309-441 lbs	170 / 67"
01-45600004	ThoraxSling	polyester, with seat supp.	S	30-60 / 66-132 lbs	94 / 37"
01-45600006	ThoraxSling	polyester, with seat supp.	M	50-80 / 110-176 lbs	104 / 41"
01-45600007	ThoraxSling.	polyester, with seat supp.	L	70-110 / 154-242 lbs	124 / 49"
01-45600008	ThoraxSling	polyester, with seat supp.	XL	100-150 / 220-330 lbs	146 / 57"
01-45690004	ThoraxSling	disposable, with seat supp.	S	30-60 / 66-132 lbs	94 / 37"
01-45690006	ThoraxSling	disposable, with seat supp.	M	50-80 / 110-176 lbs	104 / 41"
01-45690007	ThoraxSling.	disposable, with seat supp.	L	70-110 / 154-242 lbs	124 / 49"
01-45690008	ThoraxSling	disposable, with seat supp.	XL	100-150 / 220-330 lbs	146 / 57"

Rising/Standing with RiseUp600EE



Art. No.	Product name	Material	Size	Rec user weight (kg)	Lift vest width (cm)
01-45200004	StandUpSling	polyester	S	40-65 / 88-143 lbs	109 / 43"
01-45200006	StandUpSling	polyester	M	55-80 / 121-176 lbs	121 / 48"
01-45200007	StandUpSling	polyester	L	75-110 / 165-242 lbs	122 / 48"
01-45200008	StandUpSling	polyester	XL	100-160 / 220-353 lbs	136 / 54"
01-45200009	StandUpSling	polyester	XXL	150-220 / 331-485 lbs	166 / 65"


Standing/Walking



Art No.	Product name	Material	Size	Rec user weight (kg)	Lift vest width (cm)	Lift vest overall height (cm)
01-44503003	StandingVest	polyester	XS	15-40 / 33-88 lbs	82 / 32"	68 / 27"
01-44503004	StandingVest	polyester	S	20-50 / 44-110 lbs	96 / 38"	73 / 29"
01-44503006	StandingVest	polyester	M	40-60 / 88-132 lbs	110 / 43"	73 / 29"
01-44503007	StandingVest	polyester	L	50-80 / 110-176 lbs	122 / 48"	73 / 29"
01-44503008	StandingVest	polyester	XL	70-100 / 154-220 lbs	130 / 51"	73 / 29"



01-44200002	WalkingVest	polyester	XXS	15-40 / 33-88 lbs	60 / 24"	67 / 26"
01-44200003	WalkingVest	polyester	XS	25-50 / 55-110 lbs	70 / 28"	68 / 27"
01-44200004	WalkingVest	polyester	S	40-60 / 88-132 lbs	82 / 32"	68 / 27"
01-44200006	WalkingVest	polyester	M	55-80 / 121-176 lbs	96 / 38"	73 / 29"
01-44200007	WalkingVest	polyester	L	75-100 / 165-220 lbs	110 / 43"	73 / 29"
01-44200008	WalkingVest	polyester	XL	90-110 / 198-242 lbs	122 / 48"	73 / 29"
01-44200009	WalkingVest	polyester	XXL	100-130 / 220-287 lbs	142 / 56"	75 / 30"



01-47103010	BariVest	polyester	XXXL	180-300 / 397-661 lbs	190-280 / 75-110"	-
01-47103011	BariVest	polyester	XXXXL	250-570 / 551-1257 lbs	230-310 / 91-122"	-

PLEASE NOTE!

To determine the most appropriate model, size and material for the individual user, trial fitting of lifting slings is always essential. The size guide's measurements and recommendations are only intended as an aid to selecting lifting slings. The lifting sling's seat width gives an indication of how wide the user's hips (trochanter width) can be. The user's trochanter width must be within +/- 10 cm/3.9" of the lifting sling's seat width. The lift vest's width gives an indication of the girth of the user's trunk. The user's trunk girth must be within +/- 10 cm/3.9" of the lift vest's width. The lifting sling is too large if the user's seat slides out too far at the lower opening, or if the leg supports slide forward towards the hollows of the user's knees. The lifting sling is too small if the user's face is too close to the slingbar, if the leg supports chafe at the groin, or if the sling appears to be too tight or feels too tight.

Lifting slings

Accessories



AdjustmentLoop

AdjustmentLoop is an extension strap, mainly used to extend the sling's lift straps to the desired length; e.g., to enable a more recumbent position or a more upright seated position, or to achieve better balance, for example, for users with hemiplegia. AdjustmentLoop is available in two lengths, both for max 300 kg/660 lbs.

Art. No.	Product name	Description	Unit
01-49182006	AdjustmentLoop	18 cm/ 7.09", provides an extension of about 12 cm/ 4.72"	1 Pair
01-49182007	AdjustmentLoop	28 cm/ 11.02", provides an extension of about 22 cm/ 8.66"	1 Pair



HeadSupport

HeadSupport is an accessory for lifting slings without head support. HeadSupport is placed around the user's head and then connected to the sling bar hooks to provide extra support for the head. HeadSupport has three alternative loops on both sides to enable choice of the most appropriate length. We recommend analysis of different sling models to determine the need for a sling that provides head support. Max weight: 300 kg/660 lbs.

Art. No.	Product name	Description	Unit
01-49302000	HeadSupport	Polyester, one size	1 Pce



SoftLegSupport

Pressure-sensitive users sometimes need extra soft leg supports. SoftLegSupports are padded covers that can be applied to the leg sections of slings with divided leg supports. SoftLegSupports are sold in pairs and are easy for the caregiver to apply.

Art. No.	Product name	Description	Unit
01-49373003	SoftLegSupport	X-Small	1 Pair
01-49373006	SoftLegSupport	Small/Medium	1 Pair
01-49373007	SoftLegSupport	Large/X-Large	1 Pair



Lifting slings

Accessories



StrapsPadding

StrapsPadding are padded strap covers which can be used as extra padding on the lift straps of a lifting sling to alleviate pressure for users who experience discomfort. StrapsPadding can be placed in many different ways to increase the comfort, for example for users with a strongly leaning seated position and/or uncontrolled movements. They are fastened with velcro and are very easy to fit, move or remove for cleaning. StrapsPadding can be used as an accessory for all lifting slings with divided leg supports.

Art. No.	Product name	Description	Unit
01-49373300	StrapsPadding	Padded lift strap covers	2 Pcs



NameTag

NameTag is a textile marking label which is easily attached to, removed from and moved between products. It is safely attached and can be washed at up to 90° C/194°F.

Art. No.	Product name	Description	Unit
01-49700001	NameTag	Textile marking label	10 Pcs



Safety information for transfer aids.



Visual inspection

Ensure that the material and the seams are not damaged. Verify that the material is not worn or bleached.



Mechanical load or stress

Apply pressure to the buckles and handles and verify that they will not open or break under heavy loads.



Stability test

Carefully bend the slab to verify that it is elastic but still feels firm and solid.



Friction test

Low friction - Ensure that the surface is smooth and slides effortlessly. Compare with a new product.

High friction - Ensure that the surface is rough and prevents the transfer aid from moving. Compare with a new product.


Make sure your transfer aid is in perfect working order!

Transfer aids get older too. Conduct regular checks to ensure that they are safe and in perfect working order! An aid that does not meet the required demands must be removed from service immediately.

Aids with handles and straps




Check safety!



Aids with buckles




Check safety!



Aids with a stiff core




Check safety!



Aids with low or high friction




Check safety!



Aids with loops



Check safety!



Please visit our website to download all documentation and to read more about the function and maintenance of our products. On our website, you can also find video clips and instructions on how to use the products. But please don't forget to read the manual!