









WallRack is a wall-mounted hanger that is used mainly for storing disposable covers for EasySlide and RollerSlide, but can also be used for storing the products themselves.

SystemRoMedic is a series of different models of transfer assistive devices. The choice of assistive devices depends on the transfer situation and the patient's functional capability and needs.

Functional inspection


Visual inspection

Check the condition and function of the assistive device regularly. If there are signs of wear, the product must be discarded.

$\mathbf{\widetilde{1}}$ Always read the manual

Always read the manuals for all assistive devices used during a transfer. Keep the manual where it is accessible to users of the product.

Wallrack



WallRack can be mounted either right-hung or left-hung on the wall.

Easily disassembled for remounting hangers on either the left or right side of the mid-section.

Care of the product:

Clean with dishwashing liquid or disinfectant. Wipe with a cloth.

Medical Device Class I

The product meets the requirements of the Medical Devices Directive 93/42/EEC.

Contact your local distributor if you have any questions about the product and its use. See www.handicare.com for a complete list of distributors. Always make sure that you have the right version of the manual. The most recent editions of manuals are available for downloading from our website, www.handicare.com.

System**RoMedic**™

For 25 years SystemRoMedic have applied themselves wholeheartedly to developing smart and easy-to-use assistive devices for easy transfers and to making life and work easier for both patients and personnel in the care sector.

Experience, innovation and training are the basis for SystemRoMedic, a total solution for every imaginable transfer situation.

Transfer: products for moving patients between locations. Positioning: products for repositioning in the same location. Support: products that provide mobility support. Lifting: products adapted for lifting.

The philosophy behind SystemRoMedic is to prevent occupational injuries while improving the patient's sense of independence and dignity. Through a combination of training and a complete range of transfer-assistive devices, SystemRoMedic offers the means to improving both the work environment and the quality of care while enabling significant cost savings.

Their mission, to help people, has always been, and will continue to be, the driving force of innovation. They love easy transfers.



info@amilakesouthern.co.uk